

Product datasheet for **RG234689**

CDHH (CDH13) (NM_001220489) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDHH (CDH13) (NM_001220489) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CDHH
Synonyms:	CDHH; P105
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG234689 representing NM_001220489
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCAGCCGAGAACTCCGCTCGTTCTGTGCGTTCTCCTGTCCCAGGTGCTGCTAACATCTGCAGAAG
ATTTGGACTGCACTCCTGGATTTTCAGCAGAAAGTGTCCATATCAATCAGCCAGCTGAATTCATTGAGGA
CCAGTCAATTCTAAACTTGACCTTCAGTACTGTAAAGGAAACGACAAGCTACGCTATGAGGTCTCGAGC
CCATACTTCAAGGTGAACAGCGATGGCGGCTTAGTTGCTCTGAGAAACATAACTGCAAGTGGGCAAAACTC
TGTTCTGTCATGCACGGACCCCCATGCGGAAGATATGGCAGAACTCGTGATTGTCGGGGGAAAGACAT
CCAGGGCTCCTTGCAAGTAGTCGATAGTGACAGGCCAGAAAGTCCAAGTCCGGCTCACTGAAAAGGGA
GTGGATCAAGAGCCTAAAGGAATTTTCAGAATCAATGAGAACACAGGGAGCGTCTCCGTGACACGGACCT
TGGACAGAGAAGTAATCGCTGTTTATCAACTATTTGTGGAGACCACTGATGTCAATGGCAAAACTCTCGA
GGGGCCGGTGCCTCTGGAAGTCATTGTGATTGATCAGAATGACAACCGACCGATCTTTCCGGAAGGCCCC
TACATCGGCCACGTCATGGAAGGGTCACCCACAGGCACACAGTATGCGGATGACAGCCTTTGATGCGAG
ATGACCCAGCCACCGATAATGCCCTCCTGCGGTATAATATCCGTGACGAGACGCCTGACAAGCCATCTCC
CAACATGTTTACATCGATCCTGAGAAAGGAGACATTGTCAGTGTGTGTGACCTGCGCTGCTGGACCGA
GAGACTCTGAAAATCCCAAGTATGAACTGATCATCGAGGCTCAAGATATGGCTGGACTGGATGTTGGAT
TAACAGGCACGGCCACAGCCACGATCATGATCGATGACAAAAATGATCACTCACCAAAATTCACCAAGAA
AGAGTTTCAAGCCACAGTCGAGGAAGGAGCTGTGGGAGTTATTGTCAATTTGACAGTTGAAGATAAGGAT
GACCCACACAGGTGCATGGAGGGCTGCCTACCCATCATCAACGAAACCCCGGGCAGAGCTTTGAAA
TCCACACCAACCCTCAAACCAACGAAGGGATGCTTTCTGTTGTCAAACCATGGACTATGAAATTTCTGC
CTTCCACACCCTGCTGATCAAAGTGAAAAATGAAGACCCACTCGTACCCGACGTCCTACGGCCCCAGC
TCCACAGCCACCGTCCACATCACTGTCTGGATGTCAACGAGGGCCAGTCTTCTACCCAGACCCCATGA
TGGTGACCAGGCAGGAGGACCTCTCTGTGGCAGCGTGTCTGACAGTGAATGCCACGGACCCCGACTC
CCTGCAGCATCAAACCATCAGGTATTCTGTTTACAAGGACCCAGCAGGTTGGCTGAATATTAACCCCATC
AATGGGACTGTTGACACCACAGCTGTCTGGACCGTGAGTCCCCATTTGTCGACAACAGCGGTACACTG
CTCTCTTCTGGCAATTGACAGTGGCAACCCTCCCGCTACGGGCACTGGGACTTTGCTGATAACCCTGGA
GGACGTGAATGACAATGCCCGTTCATTTACCCACAGTAGCTGAAGTCTGTGATGATGCCAAAAACCTC
AGTGTAGTCATTTGGGAGCATCAGATAAGGATCTTACCCGAATACAGATCCTTTCAAATTTGAAATCC
ACAACAAGCTGTTCTGATAAAGTCTGGAAGATCTCAAGATCAACAATACACAGCCCTGGTAAGCCT
TCTTCAAATCTGAACAAGCAAACCTACAACCTGCCCATCATGGTGACAGATTCAGGGAAACCCCATG
ACGAATATCACAGATCTCAGGGTACAAGTGTGCTCCTGCAGGAATTCAAAGTGGACTGCAACGCGGCAG
GGCCCTGCGCTTCAGCCTGCCCTCAGTCTGCTCCTCAGCCTTTCAGCTTAGCTTGTCTG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG234689 representing NM_001220489
 Red=Cloning site Green=Tags(s)

MQPRTPLVLCVLLSQVLLLTSAEDLDCTPGFQQKVFHINQPAEFIEDQSILNLTFSDCCKGNDKLRYEVSS
 PYFKVNSDGGVALRNITAVGKTLFVHARTPHAEDMAELVIVGGKDIQGSLOVVDSDRPERSKFRLTGKG
 VDQEPKGI FRINENTGSVSVTRTLDREVI AVYQLFVETTDVNGKTLGEPVPLEVIDQNDNRPIFREGP
 YIGHVMEGSPGTGTTVMRMTAFDADDPATDNALLRYNIRQQTPDKPSPNMFYIDPEKGDIVTVVSPALLDR
 ETLENPKYELII EAQDMAGLDVGLTGTATATIMIDDKNDHSPKFTKKEFQATVEEGAVGVI VNLTVEDKD
 DPTTGAWRAAYTI INGNPGQSFEIHTNPQTNEGMLSVVKPLDYEISAFHTLLIKVENEDPLVPDVSYGPS
 STATVHITVLDVNEGPVFYDPMMVTRQEDLSVGSVLLTVNATDPDSLQHQITIRYSVYKDPAGWLNINPI
 NGTVDTTAVLDRESFVDNSVYALFLAIDSGNPPATGTGTLITLEDVNDNAPFIYPTVAEVCDDAKNL
 SVVILGASDKDLHPNTDPFKFEIHKQAVPDKVWKISKINNTHALVSLQLNLNKANYNLPIMVTDSGKPPM
 TNITDLRVQVCSRNKSKVDCNAAGALRFLPSVLLLSLFLACL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001220489

ORF Size: 2022 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001220489.2](#)

RefSeq Size: 3911 bp

RefSeq ORF: 2025 bp


Locus ID: 1012

UniProt ID: [P55290](#)

Cytogenetics: 16q23.3

Gene Summary: This gene encodes a member of the cadherin superfamily. The encoded protein is localized to the surface of the cell membrane and is anchored by a GPI moiety, rather than by a transmembrane domain. The protein lacks the cytoplasmic domain characteristic of other cadherins, and so is not thought to be a cell-cell adhesion glycoprotein. This protein acts as a negative regulator of axon growth during neural differentiation. It also protects vascular endothelial cells from apoptosis due to oxidative stress, and is associated with resistance to atherosclerosis. The gene is hypermethylated in many types of cancer. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2011]

Product images:



Circular map for RG234689