

Product datasheet for **RG234697**

FOXP1 (NM_001244814) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOXP1 (NM_001244814) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FOXP1
Synonyms:	12CC4; hFKH1B; HSPC215; MFH; QRF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)


**ORF Nucleotide
Sequence:**

>RG234697 representing NM_001244814
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGCAAGAATCTGGGACTGAGACAAAAAGTAAACGGTTCAGCCATCCAGAATGGGTGGGCGGCAGCA
 ACCACTTACTAGAGTGCGGCGGTCTTCGGGAGGGCGGTCCAACGAGAGACGCCGGCGTGGACATCGG
 GGAGCTGACCTCGCCACGCCAGCAGCAGCAACAGGCACTTCAGGTGGCAAGACAGCTCCTTCTT
 CAGCAGCAACAGCAGCAGCAAGTTAGTGGATTAATCTCCCAAGAGGAATGACAAAACAACCAGCTCTTC
 AGGTTCCCGTGTGAGTGGCTATGATGACACCTCAAGTTTCACTCCCAAGCAATGCAGCAGATCTCCA
 GCAACAAGTGCTGAGCCCTCAGCAGCTCCAGGTTCTCTCCAGCAGCAGCAGGCCCTCATGCTTCAACAG
 CAGCAGCTTCAAGAGTTTATAAAAAACAACAGGAACAGTTGCAGCTTCAACTTTTACAACAACAACATG
 CTGAAAAACAGCCTAAAGAGCAACAGCAGGTGGCTACCCAGCAGTTGGCTTTTCAGCAGCAGCTTTTACA
 GATGCAGCAGTTACAGCAGCAGCACCTCTGTCTTTCAGCGCCAAGGCCTTCTGACAATTCAGCCCGGG
 CAGCCTGCCCTTCCCTTCAACCTTGTCTCAAGGCATGATTCCAACAGAACTGCAGCAGCTCTGAAAAG
 AAGTGACAAGTGCTCATACTGCAGAAGAAACCACAGGCAACAATCACAGCAGTTTGGATCTGACCACGAC
 ATGTGTCTCCTCTCTGCACCTTCCAAGACCTCTTAATAATGAACCCACATGCCTTACCAATGGACAG
 CTCTCAGTCCACACTCCCAAAAGGGAAAGTTTGTCCCATGAGGAGCACCCCATAGCCATCTCTCTATG
 GACATGGTGTATGCAAGTGGCCAGGCTGTGAAGCAGTGTGCAAGATTTCCAATCATTTCTAAAACATCT
 CAACAGTGAGCATGCGCTGGACGATAGAAGTACAGCCCAATGTAGAGTACAAATGCAGGTTGTACAGCAG
 TTAGAGCTACAGCTTGCAAAAGACAAAGAACGCCTGCAAGCCATGATGACCCACCTGCATGTGAAGTCTA
 CAGAACCCAAAGCCGCCCTCAGCCCTTGAATCTGGTATCAAGTGTCACTCTCTCAAGTCCGCATCGGA
 GGCTTCTCCACAGAGCTTACCTCATACTCCAACGACCCCAACCGCCCCCTGACTCCCGTCACCCAAGGC
 CCTCTGTCAACAACCACAGCATGCACAGGTGGGACCCATCCGCAGGCGGTACTCAGACAAATACA
 ACGTGCCCATTTTCGTCAGCAGATATTGCGCAGAACCAAGAATTTTATAAGAACGCAGAAGTTAGACCACC
 ATTTACATATGCATCTTTAATTAGGCAGGCCATTCTCGAATCTCCAGAAAAGCAGCTAACACTAAATGAG
 ATCTATAACTGGTTCACACGAATGTTTGTCTTACTTCCGACGCAACCGGCCACGTGGAAGAATGCAGTGC
 GTCATAATCTTAGTCTTCAAGTGTGTTGTCGAGTAGAAAACGTTAAAGGGGAGTATGGACAGTGG
 TGAAGTAGAATTCAAAAACGAAGGCCAAAAAGATCAGTGGTAAACCTTCCCTTATTAACAAATGCAG
 AGCAGCCACGCCTACTGCACACCTCTCAATGCAGCTTTACAGGCTTCAATGGCTGAGAATAGTATACCTC
 TATACACTACCGTTCCATGGGAAATCCCACTCTGGGCAACTTAGCCAGCGCAATACGGGAAGAGCTGAA
 CGGGGCAATGGAGCATACCAACAGCAACGAGAGTGACAGCAGTCCAGGCAGATCTCCTATGCAAGCCGTG
 CATCTGTACACGTCAAAGAAGAGCCCCTCGATCCAGAGGAAGCTGAAGGGCCCCCTGTCTTAGTGACAA
 CAGCCAACCACAGTCCAGATTTTGACCATGACAGAGATTACGAAGATGAACCAAGTAAACGAGGACATGGA
 G

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Plasmid Map:


ACCN: NM_001244814

ORF Size: 2031 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001244814.3](#)

RefSeq Size: 6899 bp

RefSeq ORF: 2034 bp

Locus ID: 27086

UniProt ID: [Q9H334](#)

Cytogenetics: 3p13

Protein Families: Transcription Factors

Gene Summary: This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Forkhead box P1 protein contains both DNA-binding- and protein-protein binding-domains. This gene may act as a tumor suppressor as it is lost in several tumor types and maps to a chromosomal region (3p14.1) reported to contain a tumor suppressor gene(s). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]