

Product datasheet for **RG234713**

FAM13A (NM_001265578) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM13A (NM_001265578) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FAM13A
Synonyms:	ARHGAP48; FAM13A1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG234713 representing NM_001265578
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTTGTGAAATCATGCCTCTGCAAAGTGCTCATGTACCCCAAGTCAGCAATGTGTCTGCAACCGGAG
 AACTCTTAGAAAGAACCATCCGATCAGCTGTAGAACAACATCTTTTGTGTTAATAACTCTGGAGGTCA
 AAGTTCAGAGGACTCAGAATCTGGAACACTATCAGCATCTTCTGCCACATCTGCCAGACAGCGCCCGCGC
 CAGTCCAAGGAGCAGGATGAAGTTCGACATGGGAGAGACAAGGGACTTATCAACAAAGAAAATACTCCTT
 CTGGGTTCACACCCTTGATGATTGTATTTGAATACTCAGGAAGTCGAAAAGGTACACAAAAATACTTT
 TGGTTGTGCTGGAGAAAGGAGCAAGCCTAAACGTCAGAAATCCAGTACTAACTTTCTGAGCTTCATGAC
 AATCAGGACGGTCTTGGAATATGGAAGTCTCAATCCACACGATCTCATGAGAGAAGTGGACCTGATG
 ATTTTGAATGGATGTCTGATGAAAGGAAAGGAAATGAAAAGATGGTGGACACACTCAGCATTTTGAGAG
 CCCCACAATGAAGATCCAGGAGCATCCAGCCTATCTGACACCAACAGCAGAGAAATCAAGATGCCGGT
 GACCAGGAGGAGAGCTTTGTCTCGAAGTGCCCCAGTCCGACCTGACTGCATTGTGTGATGAAAAGAACT
 GGAAGAGCCTATCCCTGCTTTCTCTCTGCGAGCGGGAGAACAGTGACTCTGATGAAGCCCACCTCTC
 GCCGACGGCTGGGCGCCTGATCCGTGAGCTGTGGACGAAGACAGCGACCCCATGCTCTCTCTCGGTTT
 TACGCTTATGGGAGAGCAGGCAATACCTGGATGACACAGAAGTGCCTCCTTCCCAACAACTCCCAT
 CTTTCATGAGGCGGCGAAGCTCCTCTCTGGGTCTATGATGATGAGCAAGAGGACCTGACACCTGCCCA
 GCTCACACGAAGGATTCAGAGCCTTAAAAAGAAGATCCGAAGTTTGAAGATAGATTGGAAGAAGAGAAG
 AAGTACAGACCTTCCACAGTGACAAAGCAGCAATCCGGAGTTCTGAAATGGACAAATGACCTTGCCA
 AATTCGGAGACAACCTTAAAGAATCAAACTAAAGATATCTGAAGAGGACCTAACTCCAGGATGCCGCA
 GCGAAGCAACACACTCCCAAGAGTTTGGTTCCCACTTGAGAAAGAAGATGAGAAGAAGCAAGAGCTG
 GTGGATAAAGCAATAAAGCCAGTGTTGAAGCCACATTGGAATCTATTCAGAGGAAGCTCCAGGAGAAGC
 GAGCGAAAGCAGCGCCCTGAGGACATTAAGGATATGACCAAGACCAGATTGCTAATGAGAAAGTGGC
 TCTGCAGAAAGCTCTGTTATATTATGAAAGCATTATGAGCGGCCGTAACAAAGAAGCAACGGCAGGTG
 ATGAAGCCACTATACGACAGGTACCGCTGGTCAAACAGATCCTCTCCGAGCTAACACCATACCCATCA
 TTGGTTCCCCCTCAGCAAGCGGAGAAGCCCTTGTGTCAGCCAATTATCGAGGGCGAACTGCTTCCTT
 CTTCAAGGAGATAAAGGAAGAAGAGGAGGGTCAAGAGCATAGCAATGTGAAGCCAGACTTCATGGTC
 ACTCTGAAAACCGATTTTCAGTGCAGATGCTTTCTGGACCAATTCGAAGATGACGCTGATGGATTTATTT
 CCCCAATGGATGATAAAATACCATCAAAATGCAGCCAGGACACAGGGCTTCAAACTCTCCATGCTGCCTC
 AATACCTGAACTCCTGGAACACCTCCAGGAAATGAGAGAAGAAAAGAAAGGATTGGAAGAAATTCGG
 GATTTTGAAGACAACCTTTTCAGACAGAATGGAAGAAATGTCCAGAAGGAAGACCGCACTCCTATGGCTG
 AAGAATACAGTGAATATAAGCACATAAAGGCGAACTGAGGCTCCTGGAGGTGCTCATCAGCAAGAGAGA
 CACTGATTCCAAGTCCATG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG234713 representing NM_001265578
 Red=Cloning site Green=Tags(s)


MACEIMPLQSAHVPQVSNVSATGELLERTIRSAVEQHLFDVNNSSGGQSSSESGTLSASSATSARQRRR
 QSKEQDEVHRGRDKGLINKENTPSGFNHLDDCILNTQEVEKVHKNTFGCAGERSKPKRQKSSTKLSELHD
 NQDGLVNMESLNSTRSHERTGPDDFEWMSDERKGNEKDGGHTQHFESPTMKIQEHPSLSDTKQQRNQDAG
 DQEESFVSEVPQSDLTALCDEKNWEEPIPAFSSWQRENSDSDEAHLSPQAGRLIRQLLEDSDPMLSPRF
 YAYGQSRQYLDDETEVPPSPPNHSFMRRRSSSLGSYDDEQEDLTPAQLTRRIQSLKKKIRKFEDRFEEEEK
 KYRPSHSDKAANPEVLKWTNDLAKFRRQLKESKLKISEEDLTPRMRQRSNTLPKSFSGSQLEKEDEKKQEL
 VDKAIKPSVEATLESIQRKLQEKRAESSRPEDIKDMTKDQIANEKVALQKALLYESIHDRPVTKNERQV
 MKPLYDYRYLVKQILSRANTIPIIGSPSSKRSPLLQPIIEGETASFFKEIKEEEEGSDDSNVDPDFMV
 TLKTDFSARCFLDQFEDDADGFISPMDDKIPSKCSQDTGLSNLHAASIPELLEHLQEMREEKKRIRKKLR
 DFEDNFRQNGRNVQKEDRTPMAEEYSEYKHIAKALRLLEVLISKRDTSKSM

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001265578

ORF Size: 2049 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001265578.1, NP_001252507.1

RefSeq Size: 4988 bp


RefSeq ORF: 2052 bp

Locus ID: 10144

UniProt ID: O94988

Cytogenetics: 4q22.1

Product images:



Circular map for RG234713