

Product datasheet for **RG234764**

SEC14L1 (NM_001204410) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SEC14L1
Synonyms:	PRELID4A; SEC14L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG234764 representing NM_001204410
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGTGCAGAAATACCACTCCCACTGAGGGTGTACAATACCCCTTTGAATTAATTATGGCTGCCTATG
 AAAGGAGGTTCCCTACATGTCTTTGATTCCGATGTTCTGGGCGAGTACACTGTGAATGAATCAAGAG
 CGAAGATGGGGCTATTCATGTCAATTGAAAGGCGCTGCAAGCTGGATGTAGATGCACCCAGACTGCTGAAG
 AAGATTGCAGGAGTTGATTATGTTTATTTGTCCAGAAAACTCACTGAATTCGCGGAACGTACTTTGC
 ACATTGAGGCTTATAATGAAACGTTTCCAATCGGGTCATCATTAAAGAGCATTGCTACACCGTTCA
 CCTGAAAAATGAAGATTGGACCTGTTTGAACAGTCTGCAAGTTTAGATATTAATCTTTCTTTGGTTTT
 GAAAGTACAGTGGAAAAATTGCAATGAAACAATATACCAGCAACATTAAGAAAGGAAAGGAATCATCG
 AATACTACCTTCGCCAATTAGAAGAAGAGGCATAACCTTTGTGCCCCGTTGGAGTCCGCTTCCATCAC
 GCCCTCTTCAGAGACATCTTCATCATCTCCAAGAAACAAGCAGCGTCCATGGCCGTCGTCATCCAGAA
 GCTGCCCTCAAGGAGGGGCTGAGTGGTGTGCCCTCAGCAGCCCCAGCGACCTGAGCCCGTGGTGGGCA
 CCCCTGACGACAACTAGATGCCGACTACATCAAGAGATACCTGGGCGATTGACTCCGCTGCAGGAGAG
 CTGCCCTATTAGACTTCGCCAGTGGCTCCAGGAGACCCACAAGGGCAAAATCCAAAAGATGAGCATATT
 CTTGCGTTCCTCCGTGCACGGGATTTAATATTGACAAAGCCAGAGAGATCATGTGTGAGTCTTTGACGT
 GGAGAAAGCAGCATCAGGTAGACTACATTCTTGAACCTGGACCCCTCCTCAGGTCCTTCAGGATTACTA
 CGCGGGAGGCTGGCATCATCAGCACAAGATGGGCGGCCCTCTACGTGCTCAGGCTGGGGCAGATGGAC
 ACCAAGGGCTTGGTGAAGAGCGCTCGGGGAGGAAGCCCTGCTGAGATACGTTCTCTCCATAAATGAAGAAG
 GGCTAAGGCGATGCGAAGAGAATACAAAAGTCTTTGGTGGCCTATCAGCTCATGGACCTGCCTGGTGA
 CTTGGAAGGGCTGAACATGCGCCACTTGTGGAGACCTGGTGTGAAAGCGCTGCTGCGGATCATCGAGGTG
 GTGGAGGCCAACTACCTGAGACACTGGGCGGCTTCTCATCTGCGGGCGCCAGGGTATTTCTGTGC
 TCTGGACGCTGGTGTAGTCCGTTTATTGATGACAACACCAGAAGGAAGTTCTCTCATTTATGACGAAATGA
 TCACAGGGTCTCGAGGCTGCTGGATTACATCGACAAAGAGATTATTCAGATTTCTGAGTGGGGAG
 TGCATGTGCGAAGTGCCAGAGGGTGGACTGGTCCCAAAATCTCTGTACCGGACTGCAGAGGAGCTGGAGA
 ACGAAGACCTGAAGCTCTGGACTGAGACCATCTACAGTCTGCAAGCGTCTTCAAAGGAGCCCCACATGA
 GATTCTCATTCAGATTGTGGATGCCTCGTCAGTCATCACTTGGGATTTGACGTGTGCAAAGGGGACATT
 GTGTTTAACATCTATCACTCCAAGAGTGCACACAACCACCAAAAGGACTCCCTGGGAGCCACAGCA
 TCACCTCTCCGGGTGGGAACAATGTGAGCTCATAGACAAAGTCTGGCAGCTGGGCCGCGACTACAGCAT
 GGTGGAGTCGCTCTGATCTGCAAAGAAGGAGAAAGCGTGCAGGGTTCCATGTGACAGGTGGCCGGGC
 TTCTACATCTGCACTGGAATTCACAGCATGCCTGCGTGCGCCGCCAGCAGCTTCCCGGGTGGAGC
 ACGTGCTTGCCTCCCTGCAGGTCTCTTCGACAAAGTGTAAAGTGTACTACACCGAGGTGATCGGCTC
 GGAGGATTTGAGAGTTCCATGACGAGCTGGAGTCCAGCCACAGCGGCTTCTCCAGCTGAGTGCCGCC
 ACCACCTCTCCAGCCAGTCCCACTCCAGCTCCATGATCTCCAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG234764 representing NM_001204410
 Red=Cloning site Green=Tags(s)

MVQKYQSPVRVYKYPFELIMAAYERRFPTCLIPMFVGS DTVNEFKS EDGAIHVIERRCKLDVDAPRLK
 KIAGVDYVYFVQKNSLNSRERTLHIEAYNETFSNRVINEHCCYTVHPENEDWTCFEQSASLDIKSFFGF
 ESTVEKIAMKQYTSNIKKGKEIIEYYLRQLLEEIGITFVPRWSPSITPSSETSSSSSKKQAASMAVVIPE
 AALKEGLSGDALSSPAPEPVVGTDDKLDADYIKRYLGDLTPLQESCLIRLRQLQETHKGKIPKDEHI
 LRFLRARDFNIDKAREIMCQSLTWRKQHQVDYILETWPPQVLQDYAGGWHHDKDGRPLVYLRGQMD
 TKGLVRALGEEALLRYVLSINEEGLRCEENTKVFGRPISSWTCLVDLEGLNMRHLWRPGVKALLRIIEV
 VEANYPETLGRLLILRAPRVFVWLTVLSPFIDNTRRKFLIYAGNDYQGPGLLDYIDKEIIPDFLSGE
 CMCEVPEGGLVPKSLYRTAELENE DLKLTETIYQSASVFKGAPHEILIQIVDASSVITWDFDVCKGDI
 VFNIYHSKRSPQPPKKDSLGAHSITSPGGNNVQLIDKVVQLGRDYSMVESPLICEGESVQGSHTVRWPG
 FYILQWKFSMPACAASSLPRVDDVLASLQVSSHCKKVMYYTEVIGSEDFRGSMTSLESSHSGFSQLSAA
 TTSSSQSHSSSMISR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001204410

ORF Size: 2145 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001204410.1](#), [NP_001191339.1](#)

RefSeq Size: 5584 bp
 RefSeq ORF: 2148 bp
 Locus ID: 6397
 UniProt ID: [Q92503](#)
 Cytogenetics: 17q25.2-q25.3

Gene Summary: The protein encoded by this gene belongs to the SEC14 cytosolic factor family. It has similarity to yeast SEC14 and to Japanese flying squid RALBP which suggests a possible role of the gene product in an intracellular transport system. Multiple alternatively spliced transcript variants have been found for this gene; some variants represent read-through transcripts that include exons from the upstream gene C17orf86. [provided by RefSeq, Feb 2011]

Product images:



Circular map for RG234764