

Product datasheet for **RG234836**

SOX5 (NM_001261415) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SOX5 (NM_001261415) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SOX5
Synonyms:	L-SOX5; L-SOX5B; L-SOX5F; LAMSHF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234836 representing NM_001261415
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGTGTATGTCTTCCAAGCGACCAGCCTCTCCGTATGGGAAGCAGATGGAGAGGTAGCCATGGTGA
 CAAGCAGACAGAAAGTGAAGAAGAGGAGAGTGACGGGCTCCACGCTTTTACCTTCCCTTGCATGTGAG
 TTTTCCCAACAAGCCTCACTCTGAGGAATTTAGCCAGTTTCTCTGCTGACGCAAGAGACTTGTGGCCAT
 AGGACTCCCACTTCTCAGCACAATAACAATGGAAGTTGATGGCAATAAAGTTATGTCTTATTGCCCCAC
 ACAACTCATCTACCTCACCTCAGAAGGCAGAAGAAGGTGGGCGACAGAGTGGCGAGTCTTGTCTAGTAC
 AGCCCTGGGAACCTCTGAACGGCGCAAGGGCAGTTTAGCTGATGTTGTTGACACCTTGAAGCAGAGGAAA
 ATGGAAGAGCTCATAAAAACGAGCCGGAAGAAACCCCAAGTATTGAAAACTACTCTCAAAGGACTGGA
 AAGACAAGCTTCTTGAATGGGATCGGGAACTTTGGCGAAATAAAGGGACTCCCGAGAGCTTAGCTGA
 GAAAGAAAGGCAACTCATGGGTATGATCAACCAGCTGACCAGCTCCGAGAGCAGCTGTTGGCTGCCAC
 GATGAGCAGAAGAACTAGCTGCCTCTCAGATTGAGAAACAGCGTCAGCAATGGAGCTGGCCAAGCAGC
 AACAAAGACAAATGCAAGACAGCAGCAGCAGCTTCTACAGCAACAACACAAATCAATTTGCTCCAGCA
 ACAGATCCAGGTTCAAGGTGAGCTGCCGCCATTAATGATCCCGTATTCCCTCCTGATCAACGGACTG
 GCTGCAGCTGCCAGCAAGGATTCCTCCTCCCTCCAGGCTTCAGCTATAAGGCTGGATGTAGTGACCCTT
 ACCCTGTTGAGCTGATCCCACTACCATGGCAGCTGCTGCCGACGCAACACCAGGCTTAGGCCACTCCA
 ACTGCAGCAGTTATAGCTGCCAGCTAGCTGCAATGCAGGTATCTCCAGGAGGGAAGCTGCCAGGCATA
 CCCCAGGCAACCTTGGTGTCTGTATCTCCTACCAGCATTACACAGACAAGAGCACAAACAGCCCA
 CACCCAAAAGCAAGGATGAAGTGGCACAGCCACTGAACCTATCAGCTAAACCCAAGACCTCTGATGGCAA
 ATCACCCACATCACCCACCTCTCCCATATGCCAGCTCTGAGAATAAACAGTGGGGCAGGCCCTCAAA
 GCCTCTGTCCAGCAGCGTTAGCTAGTCTTACGCCAGAGTTAGCACAATAGGTTACTTAAATGACCATG
 ATGCTGTACCAAGGCAATCCAAGAAGCTCGGCAATGAAGGAGCAACTCCGACGGGAACAACAGGTGCT
 TGATGGGAAGGTGGCTGTTGTGAATAGTCTGGGTCTCAATAACTGCCAACAGAAAAGGAAAAACAACA
 CTGGAGAGTCTGACTCAGCAACTGGCAGTTAAACAGAATGAAGAAGGAAAATTTAGCCATGCAATGATGG
 ATTTCAATCTGAGTGGAGATTCTGATGGAAGTGTGGAGTCTCAGAGTCAAGAATTTATAGGGAATCCCG
 AGGGCGTGGTAGCAATGAACCCACATAAAGCGTCCAATGAATGCCTTCATGGTGTGGGCTAAAGATGAA
 CGGAGAAAAGATCCTTCAAGCCTTCTCCTGACATGCACAACCTCAACATCAGCAAGATATTGGGATCTCGT
 GGAAAGCTATGACAAACCTAGAGAAACAGCCATATTATGAGGAGCAAGCCGTCTCAGCAAGCAGCACCT
 GGAGAAGTACCCTGACTATAAGTACAAGCCAGGCCAAAGCGCACCTGCCTGGTGGATGGAAAAAGCTG
 CGCATTGGTGAATACAAGGCAATCATGCACAACAGGCGGCAGGAAATGCGGCAGTACTTCAATGTTGGGC
 AACAAAGCAGAGATCCCCATTGCCACTGCTGGTGTGTGTACCCTGGAGCCATCGCCATGGCTGGGATGCC
 CTCCCCTCACCTGCCCTCGGAGCACTCAAGCGTGTCTAGCAGCCAGAGCCTGGGATGCCTGTTATCCAG
 AGCACTTACGGTGTGAAAGGAGAGGAGCCACATATCAAAGAAGAGATACAGGCCGAGGACATCAATGGAG
 AAATTTATGATGAGTACGACGAGGAAGAGGATGATCCAGATGTAGATTATGGGAGTGACAGTGAAAAACA
 TATTGCAGGACAAGCCAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG234836 representing NM_001261415
 Red=Cloning site Green=Tags(s)


```
MSVMSSKRSPASPYGEADGEVAMVTSRQKVEEESDGLPAFHLPLHVSFPNKPHEEFQPVSLLTQETCGH
RTPTSQHNTMEVDGNKVMSSFAPHNSSTSPQKAEEGGRQSGESL SSTALGTPERRKGS LADVVDTLKQRK
MEELIKNEPEETPSIEKLLSKDWKDKLLAMGSGNFGEIKGTPESLAEKERQLMGMINQL TSLREQLLAH
DEQKLLAASQIEKQRQQMELAKQQEQIARQQQLLQQQHKINLLQQQIQVQGQLPPLMIPVFPDQRTL
AAAAQQGFLLPPGF SYKAGCSDPYVQLIPTTMAAAAAATPGLGPLQLQLYAAQLAAMQVSPGGKLPGI
PQGNLGA AVSPTS IHTDKSTNSPPPKSKDEVAQPLNLSAKPKTSDGKSPTSP TSPHMPALRINSAGPLK
ASVPAALASPARVSTIGYLNHDHAVTKAIQEARQMKELRREQVLDGKVAVVNSLGLNNCRTEKEKTT
LESLTQQLAVKQNEEGKFSHAMMDFNLSGSDSGSAGVSESR IYRESRGRGSNEPHIKRPMNAFMVWAKDE
RRKILQAFPMHNSNISKILGSRWKAMTNLEKQPYEEQARLSKQHLEKYPDYKYKPRPKRTCLVDGKKL
RIGEYKAIMRNRQEMRQYFNVGQQAQIPIATAGVYVPGA IAMAGMPSPLPSEHSSVSSSPEPGMPV IQ
STYGVKGEEPHIKEEIQ AEDINGE IYDEYDEEEDDPD VDYGSDSENHIAGQAN
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_001261415

ORF Size: 2259 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001261415.2](#)

RefSeq Size: 4357 bp

RefSeq ORF: 2262 bp

Locus ID: 6660


UniProt ID: [P35711](#)

Cytogenetics: 12p12.1

Protein Families: Transcription Factors

Gene Summary: This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The encoded protein may play a role in chondrogenesis. A pseudogene of this gene is located on chromosome 8. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG234836