

Product datasheet for **RG234837**

BRD2 (NM_001199456) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRD2 (NM_001199456) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	BRD2
Synonyms:	BRD2-IT1; D6S113E; FSH; FSRG1; NAT; O27.1.1; RING3; RNF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234837 representing NM_001199456
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTTCGGTGCCTGCTTTGCAACTTACCCTGCCAACCACCACCCCGGAGGTGTCCAATCCCAAAA
 AGCCAGGACGAGTTACCAACCAGCTGCAATACCTACACAAGGTAGTGATGAAGGCTCTGTGGAAACATCA
 GTTCGCATGGCCATTCGGCAGCCTGTGGATGCTGCAAACTGGGTCTACCGGATTATCACAAAATTATA
 AAACAGCCTATGGACATGGTACTATTAAGAGGAGACTTGAAAACAATTATTATTGGGCTGCTTCAGAGT
 GTATGCAAGATTTAATACCATGTTACCAACTGTTACATTTACAACAAGCCCACTGATGATATTGTCT
 AATGCCACAAACGCTGAAAAGATATTCCTACAGAAGTTGCATCAATGCCACAAGAAGAACAAGAGCTG
 GTAGTGACCATCCCTAAGAACAGCCACAAGAAGGGGGCCAAGTTGGCAGCGCTCCAGGGCAGTGTACCA
 GTGCCCATCAGGTGCCTGCCGTCTTTCTGTGTACACACAGCCCTGTATACTCTCCACCTGAGATACC
 TACCCTGTCTCAACATTCCCCACCCATCAGTCAATTTCTCTCCACTTCTCAAGTCTTGCCTGCTGCT
 GGACCCCGCTCCTTGCTGTTACTGCAGCTCCTCCAGCCAGCCCTTGCCAAGAAAAAGGCGTAAAGC
 GGAAAGCAGATACTACCACCCCTACACCTACAGCCATCTTGGCTCCTGGTCTCCAGCTAGCCCTCCTGG
 GAGTCTTGAGCCTAAGGCAGCACGGCTTCCCCCTATGCGTAGAGAGAGTGGTCGCCCCATCAAGCCCCCA
 CGCAAAGACTTGCTGACTCTCAGCAACAACACCAGAGCTCTAAGAAAGGAAAGCTTTCAGAACAGTTAA
 AACATTGCAATGGCATTTTGAAGGAGTACTCTCTAAGAAGCATGCTGCCTATGCTTGGCCTTTCTATAA
 ACCAGTGGATGCTTCTGCACTTGGCCTGCATGACTACCATGACATCATTAAAGACCCCATGGACCTCAGC
 ACTGTCAAGCGGAAGATGGAGAACCCTGATTACCGGGATGCACAGGAGTTTGTGCTGATGTACGGCTTA
 TGTCTCCAAGTCTATAAGTACAATCCCCAGATCACGATGTTGTGGCAATGGCACGAAAGCTACAGGA
 TGTATTTGAGTTCGTTATGCCAAGATGCCAGATGAACCACTAGAACCAGGGCCTTTACCAGTCTCTACT
 GCCATGCCCCCTGGCTTGGCCAATCGTCTTACAGTCTCCAGTGAGGAAAGTAGCAGTGAGAGCTCCT
 CTGAGGAAGAGGAGGAGGAAGATGAGGAGGACGAGGAGGAAGAAGAGAGTGAAAGCTCAGACTCAGAGGA
 AGAAAGGGCTCATCGCTTAGCAGAACTACAGGAACAGCTTCGGGCAGTACATGAACAACCTGGCTGCTCTG
 TCCCAGGGTCCAATATCCAAGCCCAAGAGGAAAAGAGAGAAAAAGAGAAAAAGAAACGGAAGGCAG
 AGAAGCATCGAGGCCGAGCTGGGCGGATGAAGATGACAAGGGGCTAGGGCACCCCGCCACCTCAACC
 TAAGAAGTCCAAGAAAGCAAGTGGCAGTGGGGTGGCAGTGTGCTTTAGGCCCTTCTGGCTTTGGACCT
 TCTGGAGGAAGTGGCACCAAGCTCCCCAAAAGGCCACAAGACAGCCCACTGCCCTGCCTACAGGTT
 ATGATTCAGAGGAGGAGGAAGAGAGCAGGCCATGAGTTACGATGAGAAGCGGCAGCTGAGCCTGGACAT
 CAACAAATTACCTGGGGAGAAGCTGGGCGGAGTTGTGCATATAATCCAAGCCAGGGAGCCCTCTTTACGT
 GATTCAAACCCAGAAGAGATTGAGATTGATTTTGAACACTCAAGCCATCCACACTTAGAGAGCTTGAGC
 GCTATGTCCTTTCTGCCTACGTAAGAAACCCCGAAGCCCTACACCATTAAGAAGCCTGTGGGAAAGAC
 AAAGGAGGAAGTGGCTTTGGAGAAAAAGCGGAATTAGAAAAGCGGTTACAAGATGTGAGCGGACAGCTC
 AATTCTACTAAAAAGCCCCCAAGAAAGCGAATGAGAAAAAGAGTATCCTCTGCACAGCAAGTAGCAG
 TGTCAGCCTTAGCGCTCCAGCTCCAGCTCAGATTCCAGCTCCTCTCTCTCTCGTCGTCGCTTTCAGA
 CACCAGTGATTCAGACTCAGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG234837 representing NM_001199456
 Red=Cloning site Green=Tags(s)


```
MASVPALQLTPANPPPPEVSNPKKGRVTNQLQYLHKVVMKALWKHQFAWPFROPVDAVKLGLPDYHKII
KQPMDMGTIKRRLNNYWAASECMQDFNTMFTNCYIYNKPTDDIVLMAQTLEKIFLQKVASMPQEEQEL
VVTIPKNSHKKGAKLAALQGSVTS AHQVPAVSSVSHTAL YTPPEIPTTVLNIPHPSVVISSPLLKSLHSA
GPPLLA VTAAPPAQPLAKKKGVKPKADTTPTPTAILAPGSPASP PGSLEPKAARLPPMRRESGRPIKPP
RKDL PDSQQQHQS S K K G L S E Q L K H C N G I L K E L L S K K H A A Y A W P F Y K P V D A S A L G L H D Y H D I I K H P M D L S
TVKRKMRNDYRDAQEF AADVRLMFSNCKYKNPPDHDVVMARKLQDVFEFRYAKMPDEPLEPGPLPVST
AMP PGLAKSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSSDSEEERHLAELQEQLRVHEQLAAL
SQGPISKPKRREKKEKKKKRKA EKHGRAGADEDKGRAPRPPQPKS K K A S G S G G S A A L G P S G F G P
S G G S G T K L P K K A T K T A P P A L P T G Y D S E E E E S R P M S Y D E K R Q L S L D I N K L P G E K L G R V V H I I Q A R E P S L R
D S N P E E I E I D F E T L K P S T L R E L E R Y V L S C L R K K P R K P Y T I K K P V G K T K E E L A L E K K R E L K R L Q D V S G Q L
N S T K K P P K A N E K T E S S S A Q Q V A V S R L S A S S S S D S S S S S S S S S S D T S D S D S G
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001199456

ORF Size: 2262 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199456.1](#), [NP_001186385.1](#)

RefSeq Size: 3210 bp

RefSeq ORF: 2265 bp

Locus ID: 6046

UniProt ID: [P25440](#)

Cytogenetics: 6p21.32

Protein Families: Protein Kinase

Gene Summary: This gene encodes a transcriptional regulator that belongs to the BET (bromodomains and extra terminal domain) family of proteins. This protein associates with transcription complexes and with acetylated chromatin during mitosis, and it selectively binds to the acetylated lysine-12 residue of histone H4 via its two bromodomains. The gene maps to the major histocompatibility complex (MHC) class II region on chromosome 6p21.3, but sequence comparison suggests that the protein is not involved in the immune response. This gene has been implicated in juvenile myoclonic epilepsy, a common form of epilepsy that becomes apparent in adolescence. Multiple alternatively spliced variants have been described for this gene. [provided by RefSeq, Dec 2010]