

Product datasheet for **RG234867**

Metabotropic Glutamate Receptor 4 (GRM4) (NM_001256814) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Metabotropic Glutamate Receptor 4 (GRM4) (NM_001256814) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Metabotropic Glutamate Receptor 4
Synonyms:	GPRC1D; mGlu4; MGLUR4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234867 representing NM_001256814
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTGCAAACCTCTTCCCAAGCTATTTCCCATGATGGCGCCAAGAGAAAAAGAGGACCCTTAGGACTT
 CAGGGCCATGCTTTGGGGTGGAGGACAGATACCCAGATCAGCTACGCCCTCCACAGCGCCAGACCTGAG
 TGACAACAGCCGCTACGACTTCTTCTCCCGCTGGTGCCCTCGGACACGTACCAGGCCAGGCCATGGTG
 GACATCGTCCGTGCCCTCAAGTGAAGTATGTGTCCACAGTGGCCTCGGAGGGCAGCTATGGTGAGAGCG
 GTGTGGAGGCCTTCATCCAGAAGTCCCCTGAGGACGGGGCGTGTGCATCGCCAGTCCGTGAAGATACC
 ACGGGAGCCCAAGGCAGGCGAGTTCGACAAGATCATCCGCCCTCCTGGAGACTTCGAACGCCAGGGCA
 GTCATCATCTTTGCCAACGAGGATGACATCAGGCGTGTGCTGGAGGCAGCACGAAGGGCCAACCAGACAG
 GCCATTTCTTCTGGATGGGCTCTGACAGCTGGGCTCCAAGATTGCACCTGTGCTGCACCTGGAGGAGGT
 GGCTGAGGGTGTGTACGATCCTCCCAAGAGGATGTCCGTACGAGGCTTCGACCGCTACTTCTCCAGC
 CGCAGCTGGACAACAACCGCGCAACATCTGGTTTCCGAGTTCTGGGAGGACAACCTCCACTGCAAGC
 TGAGCCGCCACGCCCTCAAGAAGGGCAGCCACGTCAAGAAGTGCACCAACCGTGAGCGAATTGGGCAAGG
 TTCAGCTTATGAGCAGGAGGGGAAGGTGCAAGTTTGTGATCGATGCCGTGTACGCCATGGGCCACGCGCTG
 CACGCCATGCACCGTGACCTGTGTCCCGCCGCGTGGGGCTCTGCCCGCGCATGGACCCTGTAGATGGCA
 CCCAGCTGCTTAAGTACATCCGAAACGTCAACTTCTCAGGCATCGCAGGGAACCCCTGTGACCTTCAATGA
 GAATGGAGATGCGCCTGGGCGCTATGACATCTACCAATACCAGCTGCGCAACGATTTCTGCCGAGTACAAG
 GTCATTGGCTCTGGACTGACCACCTGCACCTTAGAATAGAGCGGATGCACTGGCCGGGAGCGGGCAGC
 AGCTGCCCGCTCCATCTGCAGCTGCCCTGCCAACCAGGTTGAGCGGAAGAAGACAGTGAAGGGCATGCC
 TTGCTGTGGCACTGCGAGCCTTGACACAGGGTACCAGTACCAGGTGGACCGCTACACCTGTAAAGCGTGT
 CCCTATGACATGCGGCCACAGAGAACCACAGGGCTGCCGCCCATCCCCATCATCAAGCTTGAGTGGG
 GCTCGCCCTGGGCGTGTGCCCTCTTCTGGCCGTGGTGGGCATCGCTGCCACGTTGTTCTGGTGTGAT
 CACCTTTGTGCGTACAACGACACGCCCATCGTCAAGGCCTCGGGCCGTGAAGTACGCTACGTGCTGCTG
 GCAGGCATCTTCTGTGCTATGCCACCACCTTCTCATGATCGCTGAGCCCGACCTTGGCACCTGCTCGC
 TGCGCCGAATCTTCTGGGACTAGGGATGAGCATCAGCTATGCAGCCCTGCTACCAAGACCAACCCGAT
 CTACCGCATCTTCGAGCAGGGCAAGCGCTCGGTACGTGCCACGCTTCATCAGCCCCGCTCACAGCTG
 GCCATCACCTTACGCTCATCTCGCTGCAGCTGTGGGCATCTGTGTGTTTGTGGTGGACCCCTCCC
 ACTCGGTGGTGGACTTCCAGGACCAGCGGACACTCGACCCCGCTTCGCCAGGGTGTGCTCAAGTGTGA
 CATCTCGGACCTGTGCTCATCTGCCTGCTGGGCTACAGCATGCTGCTCATGGTACGTCACCGTGTAT
 GCCATCAAGACACGCGCGTGGCCGAGACCTTCAATGAGGCCAAGCCATTGGCTTACCATGTACACCA
 CTTGATCGTCTGGTGGCTTCCATCCCCATCTTCTTTGGCACCTCGCAGTCGGCCGACAAGCTGTACAT
 CCAGACGACGACGCTGACGGTCTCGGTGAGTCTGAGCGCCTCGGTGTCCCTGGGAATGCTCTACATGCC
 AAAGTCTACATCATCTTCCACCCGGAGCAGAACGTGCCAAGCGCAAGCGCAGCCTCAAAGCCGTCG
 TTACGGCGCCACCATGTCCAACAAGTTCACGCAGAAGGGCAACTCCGGCCCAACGGAGAGGCCAAGTC
 TGAGCTCTGCAGAACCTTGAGGCCACGCGCTGGCCACCAACAGACTTACGTCACTTACACCAACCAT
 GCAATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG234867 representing NM_001256814
Red=Cloning site Green=Tags(s)

```
MVQTLPKLFPHDGAKRKRRLRTSGPCFGGGGQIPQISYASTAPDLSDNSRYDFFSRVPSDTYQAQAMV
DIVRALKWNVYSTVASEGSYGESGVEAFIQKSREDDGGVCIQSVKIPREPKAGEFDKIIRLLETSNARA
VIIFANEDDIRRVLEAARRANQTGHFFWMGSDSWGSKIAPVLHLEEVAEGAVTILPKRMSVRGFDRYFSS
RTLDDNRRNIWF AEFWEDNFHCKLSRHALKKGSHVKKCTNRERIGQDSAYEQEGKVQFVIDAVYAMGHAL
HAMHRDLCPGRVGLCPRMDPVDGTQLLKYIRNVNFSGIAGNPVTFNENGDAPGRYDIYQYQLRNDSAEYK
VIGSWTDHLHLRIERMHWPGSGQQLPRISCSLPCQPGERKKTVKGMPCCWHCEPCTGYQYQVDRYTCKTC
PYDMRPTENRTGCRPIPIIKLEWGSWAVLPLFLAVVGI AATLFVVITFVRYNDTPIVKASGRELSYVLL
AGIFLCYATTFLMIAEPDLGTCSLRRIFLGLGMSISYAALLTKTNRIYRIFEQGRSVSAPRFISPASQL
AITFSLISLQLLGI CVWFVVDPSHSVDFDQDRTLDPRFARGVLKCDISDLSLICLLGYSMLLMVTCTVY
AIKTRGVPETFNEAKPIGFMYTTCIVWLAFIPIFFGTSQSADKLYIQTTTLTVSVLSASVSLGMLYMP
KVYIILFHPEQNVPKRKRSLKAVVTAATMSNKFTQKGNFRPNGEAKSEL CENLEAPALATKQTYVTYTNH
AI
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001256814


ORF Size: 2316 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256814.1, NP_001243743.1</u>
RefSeq Size:	3331 bp
RefSeq ORF:	2318 bp
Locus ID:	2914
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction, Taste transduction
Gene Summary:	L-glutamate is the major excitatory neurotransmitter in the central nervous system and activates both ionotropic and metabotropic glutamate receptors. Glutamatergic neurotransmission is involved in most aspects of normal brain function and can be perturbed in many neuropathologic conditions. The metabotropic glutamate receptors are a family of G protein-coupled receptors, that have been divided into 3 groups on the basis of sequence homology, putative signal transduction mechanisms, and pharmacologic properties. Group I includes GRM1 and GRM5 and these receptors have been shown to activate phospholipase C. Group II includes GRM2 and GRM3 while Group III includes GRM4, GRM6, GRM7 and GRM8. Group II and III receptors are linked to the inhibition of the cyclic AMP cascade but differ in their agonist selectivities. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]

Product images:



Circular map for RG234867