

Product datasheet for **RG234882**

CLCN3 (NM_001243374) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLCN3 (NM_001243374) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CLCN3
Synonyms:	CIC-3; CLC3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG234882 representing NM_001243374
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
 GCC

ATGGATGCTTCTCCGACCTTATTTGCCCTTATGACGGGGGAGGAGACAATATCCCTGAGGGAATTAC
 ATAAAAGAGGAACTCATTATACAATGACAAATGGAGGCAGCATTAAACAGTTCTACACATTTACTGGATCT
 TTTGGATGAACCAATTCAGGTGTTGGTACATATGATGATTTCCATACTATTGATTGGGTGCGAGAAAA
 TGTAAAGACAGAGAAAGGCATAGACGGTCAACAGCAAAAAGAAAGAAATCAGCATGGGAAATGACAAAA
 GTTTGTATGATGCGTGGTCCAGGATGGCTAGTAGTAACACTAACAGGATTGGCATCAGGGGCACTGGCCGG
 ATTAATAGACATTGCTGCCGATTGGATGACTGACCTAAAGGAGGGCATTGCTTGTAGTGGTGTGGTAC
 AACACGAAACAGTGTCTGTTGGGGATCTAATGAAACAACATTTGAAGAGAGGGATAAATGTCCACAGTGA
 AAACATGGGCAGAATTAATCATAGGTCAAGCAGAGGGTCTGGTCTTATATCATGAACTACATAATGTA
 CATCTTCTGGGCCTTGAGTTTTGCCCTTCTGTCAGTTTCCCTGGTAAAGGATTTGCTCCATATGCCTGT
 GGCCTCTGGAATCCAGAGATTAATACTATTTAAGTGGATTCATCATCAGAGTTACTTGGGAAAAATGGA
 CTTAATGATTAACCATCACATTAGTCTGGCTGTGGCATCAGGTTTGTAGTTAGGAAAAAGAGTCC
 CCTGGTACATGTTGCCTGTTGCTGCGGAAATATCTTTTCTACCTCTTTCCAAAGTATAGCACAACGAA
 GCTAAAAAAGGGAGGTGCTATCAGCTGCCTCAGCTGCAGGGTTTCTGTAGCTTTTGGTGCACCAATTG
 GAGGAGTCTTTTTAGCCTGGAAGAGGTTAGCTATTATTTTCTCTCAAACTTTATGGAGATCATTTTT
 TGCTGCTTAGTGGTGCATTTGTTTTGAGGTCCATCAATCCATTTGGTAACAGCCGTCTGGTCTTTTT
 TATGTGGAGTATCATACCCATGGTACCTTTTTGAAGTGTTCCTTTTATCTTCTAGGGGATTTGGAG
 GGCTTTGGGAGCCTTTTTATTAGGGCAATATTGCTGGTGTGTCGACGCAAGTCCAGAAATTTGG
 AAAGTATCCCGTTCTGGAAGTCATTATTGTTGCAGCCATTACTGCTGTGATAGCCTTCCCTAATCCATAC
 ACTAGGCTAAACACCACTGAACTGATCAAAGAGCTTTTTACAGACTGTGGTCCCCTGGAATCCTCTTCTC
 TTTGTGACTACAGAAATGACATGAATGCCAGTAAATGTGATGACATTCCTGATCGTCCAGCAGGCAT
 TGGAGTATATTCAGTATATGGCAGTTATGCCTGGCACTCATATTTAAATCATAATGACAGTATCACT
 TTTGGCATCAAGTTCATCAGGCTTGTTCATCCCGCAGCATGGCCATTGGAGCGATCGCAGGAAGGATTG
 TGGGGATTGCGGTGGAGCAGCTTGCTACTATCACCACGACTGGTTTATCTTTAAGGAGTGGTGTGAGGT
 CGGGGCTGATTGCATTACACCTGGCCTTATGCCATGGTGGTGTGCTGCATGCTTAGGTGGTGTGACA
 AGAATGACTGTCTCCCTGGTGGTATTGTTTTGAGCTTACTGGAGGCTTGAATATATTGTTCCCTTA
 TGCTGCAGTCATGACCAGTAAATGGGTTGGAGATGCCTTTGGCAGGGAAGGCATTTATGAAGCACACAT
 CCGATTAATGGATACCCTTTCTTGATGCAAAAAGAAATTCATCATACCACCCTGGCTGCTGACGTT
 ATGAGACCTCGAAGGAATGATCCTCCCTTAGCTGTCTGACACAGGACAATATGACAGTGGATGATATAG
 AAAACATGATTAATGAAACCAGCTACAATGGATTTCTGTGATAATGTCAAAAAGAAATCTCAGAGATTAGT
 GGGATTTGCCCTCAGAAGAGACCTGACAATTGCAATAGAAAGTCCAGGAAAAAACAAGAAGGTATCGTT
 GGCAGTTCTCGGGTGTGTTTTGCACAGCACACCCCATCTCTTCCAGCAGAAAGTCTCGGCCATTGAAGC
 TTCGAAGCATTCTGACATGAGCCCTTTACAGTGACAGACCACCCCAATGGAGATCGTGGTGGATAT
 TTTCCGAAAGCTGGGACTGAGGCAGTGCCTTGTAACTACAATGGGCGCCTCCTTGGCATTATAACAAAA
 AAAGATATCCTCCGGCATATGGCCAGACGGCAACCAAGACCCCGCTTCAATAATGTTCAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG234882 representing NM_001243374
 Red=Cloning site Green=Tags(s)

```
MDASSDPYLPYDGGGDNIPRELHKGTHYMTNGGSINSSTHLLDLLDEPIPGVGTYYDFHTIDWVREK
CKDRERHRRINSKKKESAWEMTKSLYDAWSGLVVTLTGLASGALAGLIDIAADWMTDLKEGICLSALWY
NHEQCCWGSNETTFEERDKCPQWKTWAEIIIGQAEGPGSYIMNYIMYIFWALSFAFLAVSLVKVFPYAC
GSGIPEIKTILSGFIIRGYLGKWTLMIKTITLVLAVASGLSLGKEGPLVHVACCCGNIFSYLFPKYSTNE
AKKREVLSAASAAGVSVAFGAPIGGVLFSLLEEVSYFPLKTLWRSFFAALVAAFVLRSLNPFNSRLVLF
YVEYHTPWYLFELFPFILLGVFGGLWGAFFIRANIAWCRRRKSTKFGKYPVLEVIIVAAITAVIAFPNPY
TRLNLTSELIKELFTDCGPLESSSLCDYRNDMNASKIVDDIPDRPAGIGVYSAIWQLCLALIFKIIMTVFT
FGIKVPSGLFIPMAIGAIAGRIVGIAVEQLAYYHHDWFIFKEWCEVGADCITPGLYAMVAAAACLGVT
RMTVSLVVIVFELTGGLEYIVPLMAAVMTSKWVGDAFGREGIYEAHIRLNGYPFLDAKEEFTHHTLAADV
MRPRRNDPPLAVLTQDNMTVDDIENMINETS YNGFPVIMSKESQRLVGFALRRDLTIAIESARKKQEGIV
GSSRVCF AQHTPSLPAESPRPLKLR SILDMPFTVTDHTPMEIVVDIFRKLGLRQCLVTHNGRLLGIITK
KDILRHMAQTANQDPASIMFN
```

TRTRPLE - GFP Tag - V

Restriction Sites: AscI-MluI

Cloning Scheme:



ACCN: NM_001243374


ORF Size: 2373 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001243374.1</u> , <u>NP_001230303.1</u>
RefSeq Size:	5891 bp
RefSeq ORF:	2376 bp
Locus ID:	1182
UniProt ID:	<u>P51790</u>
Cytogenetics:	4q33
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
Gene Summary:	<p>This gene encodes a member of the voltage-gated chloride channel (ClC) family. The encoded protein is present in all cell types and localized in plasma membranes and in intracellular vesicles. It is a multi-pass membrane protein which contains a ClC domain and two additional C-terminal CBS (cystathionine beta-synthase) domains. The ClC domain catalyzes the selective flow of Cl⁻ ions across cell membranes, and the CBS domain may have a regulatory function. This protein plays a role in both acidification and transmitter loading of GABAergic synaptic vesicles, and in smooth muscle cell activation and neointima formation. This protein is required for lysophosphatidic acid (LPA)-activated Cl⁻ current activity and fibroblast-to-myofibroblast differentiation. The protein activity is regulated by Ca²⁺/calmodulin-dependent protein kinase II (CaMKII) in glioma cells. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]</p>

Product images:



Circular map for RG234882