

## Product datasheet for **RG234925**

### AP180 (SNAP91) (NM\_001256717) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	AP180 (SNAP91) (NM_001256717) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AP180
Synonyms:	AP180; CALM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>RG234925 representing NM\_001256717  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAGGGCGTTCAAGAATAATAACAATTTGATCCAGGCTACCAACGAGACCAATGTTAATATTCCTCAGA  
 TGGCCGACACTCTCTTTGAGCGGCAACAACAGTAGCTGGTGGTTGTGTTAAGGCTTTAGTGACAAC  
 ACATCATCTCATGGTGCATGAAATGAGAGATTTATTCAATATTTGGCTTCTAGAATACTACTATTCAAT  
 CTCAGCAATTTTTGGACAAAAGTGATCCCATGGTTATGATATGTCTACCTTCATAAGGCGCTATAGTA  
 GATATTTGAATGAAAAGGCTTTTTCTTACAGACAGATGGCCTTTGATTTTGCCAGGGTGAAGAAAGGGC  
 CGATGGTGAATGAGGACAATGGCTCCCGAAAAGCTGCTAAAGAGTATGCCAATACTACAGGACAATTT  
 GATGCACTGCTTGAATTTGATGTGCATCCAAATGAACTAACAATGGTGCATAAATGCAGCATTTATGC  
 TTCTTTCAAAGATCTTATCAAATTTTTGCTTGTACAATGATGGTGTATTAACCTACTCGAAAAGTT  
 TTTTGAATGAAGAAAGGACAATGTAAGATGCTCTAGAAATTTACAAACGATTTCTAAGTGAATGACA  
 CGAGTGTCTGAATTTCTCAAGTTGCAGAGCAAGTTGGTATTGATAAAGGTGACATTCCTGACCTCACAC  
 AGGCTCCCAGCAGTCTTATGGAGACGCTTGAACAGCATCTAAATACATTAGAAGGAAAGAAACCTGGAAA  
 CAAATCTGGTGCCTCCCTCTCATTAAAGTAAGTCTTCTCCAGCCACAACCTGTTACGCTCCTAATTCTACA  
 CCAGCTAAAACCTATTGACACATCCCCACCGTTGATTTATTTGCAACTGCATCTGCGGCTGTCCCAGTCA  
 GCACTTCTAAACCATCTAGTGATCTCCTGGACCTCCAGCCAGACTTTTCTCTGGAGGGGAGCAGCAGC  
 CGCAGCACCAGCACCACCACCTGCTGGAGGAGCCACTGCATGGGAGACCTTTTGGGAGAGGATTCT  
 TTGGCTGCACTTTCTCTGTCCCTCTGAAGCAGAGATTCAGATCCATTTGCACCAGAACCTACCCCTC  
 TACTACAACCTGCTGAAATTTGCAACTGCCTCAGCTTCTGCCTCCACTACTACAACCTGTTACTGCTGC  
 TGCTGAAGTGGATCTCTTTGGAGATGCCTTTGCAGCTTCTCCTGGGAGGCCCTGCAGCATCCGAAGGG  
 GCCGCCGACCCAGCTACCCCAACCCCTGTAGCAGCAGCACTTGATGCATGTTTCAGGAAATGACCCCTTTG  
 CCCCCTGTAAGGTAGTGACAGAGGCTGCACCTGAGCTGGACCTCTTTGCAATGAAGCCACCTGAGACCAG  
 TGTTCTGTAGTTACCCCTACAGCTAGCACAGCCCTCCAGTTCCCGCAACTGCTCCTTCTCCTGCTCCT  
 GCCGTTGACGCTGCTGCTGCCACTACTGCTGCCACCGCCGCTGCCACCACCTACCACCACCTCCG  
 CTGCCACCGCCACCACTGCTCCTCCTGCTTAGATATCTTTGGTATGCTTTCTCCTCTCCACCACAAGG  
 GGCCTCTCCTGTGCTGAGAGTTCTCTCACTGCTGACCTCTTATCTGATGCATTTGGAAGTAGTGTCTCT  
 GAACCCCAACCTGCATCTCAGGCTGCTTCTAGTTCATCAGCATCGGCAGACCTACTAGCTGGATTTGGGG  
 GTTCTTTTCATGGCGCTTCCCCATCTCCAGTGACTCCAGCTCAGAATAACCTGCTACAGCCCAATTTTGA  
 GGCAGCTTTTGGGACAACGCCTTCAACTCCAGCAGCAGCTCCTTTGATCCATCAGTGTTTATGATGTTCTA  
 GGTGATCTTTTATGCAACCATGGCACCAGCTGGGAGCCTGCACCTGTCTCAATGGTACCACCCAGTC  
 CTGCAATGGCAGCCAGCAAAAGCCCTTGAAGTATCTTGATTCATCTCTTCCAGCTTAGTAGGCAATCT  
 TGGAAATTTCTGGTACCACAACAAAAAGGGAGATCTTCAAGTGAATGCTGGAGAGAAAAAGTTGACTGGT  
 GGAGCCAACCTGGCAGCCTAAAGTAGCTCCAGCAACCTGGTCAGCAGGCGTTCCACCAAGTGCACCTTTGC  
 AAGGAGCTGTACCTCAAACAGTTCAGTTCCTCCTGTTGCCGGGCCCCATCGTTGGACAACCTGGAGC  
 AGGATTTGGAATGCCTCCTGCTGGGACAGGCATGCCATGATGCCTCAGCAGCGGTGATGTTTGCACAG  
 CCCATGATGAGGCCCCCTTTGGAGCTGCGCTGTACCTGGCACGCAGCTTTCTCCAAGCCCTACACCTG  
 CCAGTCAGAGTCCCAAGAAACCTCCAGCAAGGACCCATTAGCGGATCTTAACATCAAGGATTTCTTG

**ACCGTACGCGGCCGCTCGAG** – GFP Tag – GTTTAA

**Protein Sequence:** >RG234925 representing NM\_001256717  
 Red=Cloning site Green=Tags(s)



```
MRAFKNNNNLIQATNETNVNIPQMADTLFERATNSSWVVVKALVTTHLMVHGNERFIQYLASRNTLFN
LSNFLDKSGSHGYDMSTFIRRYRYLNEKAFSYRQMAFDFAFVKKGADGVMRTMAPEKLLKSMPIQGGI
DALLEFDVHPNELTNGVINAAFMFLFKDLIKLFACYNDGVINLLEKFFEMKKGQCKDALEIYKRFTRMT
RVSEFLKVAEQVIGDKGIDIPDLTQAPSSLMETLEQHLNLTLEGKKPGNKSGAPSLSKSSPATTVTSPNST
PAKTIDTSPPVDL FATASAAVPVSTSKPSSDLLDLQPDFSSGGAAAAAPAPPPAGGATAWGDLGDS
LAALSSVPSEAQISDPFAPEPTPTTTAEIATASASASTTTTAVTAEVDFGDAFAASPGEAPAASEG
AAPATPTPVAALDACSGNDPFAPSEGSAAPELDFAMKPPETSVPVVTPTASTAPPVPATAPSPAP
AVAAAAATTAATAAATTTTTTSAATATTAPPALDIFGDAFSSPPQGASPVPESSLTADLLSDFGSSAS
EPQPASQAASSSSASADLLAGFGGSFMAPSPSPVTPAQNNLLQPNFEAAFGTTTPSTSSSSFDPSVFDGL
GDLLMPTMAPAGQPAPVSMVPPSPAMAASKALGSDLDSSLASLVGNLGISGTTTCKGDLQWNAGEKKLTG
GANWQPKVAPATWSAGVPPSAPLQGAVPPTSSVPPVAGAPSVGQPGAGFGMPPAGTGMPMPQPPVMFAQ
PMMRPPFGAAAVPGTQLSPSPTASQSPKKPPAKDPLADLNIKDFL
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**ACCN:** NM\_001256717

**ORF Size:** 2448 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001256717.1](#), [NP\\_001243646.1](#)

**RefSeq Size:** 4529 bp

**RefSeq ORF:** 2451 bp


**Locus ID:** 9892

**UniProt ID:** [O60641](#)

**Cytogenetics:** 6q14.2

**Gene Summary:** Adaptins are components of the adapter complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelia induces their assembly into 60-70 nm coats (By similarity).[UniProtKB/Swiss-Prot Function]

## Product images:



Circular map for RG234925