

Product datasheet for **RG235612**

C1orf150 (GCSAML) (NM_001281836) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C1orf150 (GCSAML) (NM_001281836) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: C1orf150
Synonyms: C1orf150
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235612 representing NM_001281836.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTACATTTGAAAGAAAACTTCAAGATCAAGATAAGAAAAGCCAAGAAGTTTCATCCACTTCTAAT
CAGGAAAACGAGAATGGCAGTGGTTCTGAAGAAGTGTGCTACACTGTCATTAATCACATCCCCATCAG
AGATCCTCCCTGAGCTCCAATGATGATGGCTATGAGAACATTGACTCCCTCACAAGGAAAGTGAGACAG
TTTAGAGAAAGGTCAGAGACAGAATATGCCCTTCTTAGGACTTCTGTTAGTAGGCCTTGTCTCCTGCACC
CATGAGCATGATTATGAAGTTGTGTTCCACAC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235612
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MTTFERKLQDQDKKSQEVSSSTSNQENENSGSSEEVCYTVINHIPHQRSSLSSNDDGYENIDSLTRKVRQ
FRERSETEYALLRTSVSRPCSCTHEHDYEVVFPH
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: SgfI-MluI



Cloning Scheme:




ACCN:	NM_001281836
ORF Size:	309 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001281836.2
RefSeq Size:	3978 bp
RefSeq ORF:	312 bp
Locus ID:	148823
UniProt ID:	Q5IQS6
Cytogenetics:	1q44
MW:	12.4 kDa

Gene Summary:

This gene encodes a protein thought to be a signaling molecule associated with germinal centers, the sites of proliferation and differentiation of mature B lymphocytes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]

Product images:



Circular map for RG235612