

Product datasheet for **RG235732**

VSTM1 (NM_001288793) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: VSTM1 (NM_001288793) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: VSTM1
Synonyms: SIRL-1; SIRL1; UNQ3033
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235732 representing NM_001288793.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGCAGAATTCCTCTCCCTGCTTTGCCTCGACACCAGAACCATCTTTGTCGCCATCTTCAGCTGC
ATCTCCATCCTTCTCCTCTTCTCAGTCTTCATCATCTACAGATGCAGCCAGCACAGTTCATCATCT
GAGGAATCCACCAAGAGAACCAGCCATTCCAAACTTCCGGAGCAGGAGGCTGCCGAGGCAGATTTATCC
AATATGGAAGGGTATCTCTCTCGACGGCAGACCCCAAGGAGTGACCTATGCTGAGCTAAGCACCAGC
GCCCTGTCTGAGGAGCTTCAGACACCACCCAGGAGCCCCAGGATCTCATGAATATCGGGCACTGAAA
GTG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235732
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MTAEFLSLLCLDTRTIFVAIFSCISILLFLSVFIYRCSQHSSSEESTKRTSHSKLPEQEAAEADLS
NMERVSLSTADPQGVTYAELSTALSAAASDTTQEPGPSHEYAALKV
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: Sgfl-Mlul



[View online »](#)


Cloning Scheme:


- ACCN:** NM_001288793
- ORF Size:** 348 bp
- OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- RefSeq:** [NM_001288793.1](#), [NP_001275722.1](#)
- RefSeq Size:** 674 bp
- RefSeq ORF:** 351 bp
- Locus ID:** 284415
- UniProt ID:** [Q6UX27](#)
- Cytogenetics:** 19q13.42
- Protein Families:** Transmembrane
- MW:** 13.1 kDa

Gene Summary:

Isoform 2: Behaves as a cytokine, promoting IL17A secretion by CD4+ T-cells, and differentiation and activation of IL17 producing helper T-cells (TH17).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG235732