

Product datasheet for **RG235782**

C11orf51 (ANAPC15) (NM_001278489) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C11orf51 (ANAPC15) (NM_001278489) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: C11orf51
Synonyms: APC15; C11orf51; HSPC020
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235782 representing NM_001278489.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGTCCACTTTGTTCCCCTCACTCTTCCCTCGTGTGACTGAGACTCTGTGGTTTAACTGGATCGACCC
TGTGTGGAAGAGACAGAGCTGCAGCAGCAGGAACAGCAGCATCAGGCCTGGCTCCAAAGCATCGCGGAG
AAAGACAACAACCTGGTTCCTATTGGCAAGCCAGCCTCAGAGCACTATGATGACGAGGAAGAAGAGGAT
GATGAAGATGATGAGGATAGTGAAGAGGACTCAGAGGATGATGAGGATATGCAGGACATGGACGAGATG
AATGACTACAATGAGTACCGGATGATGGAGAGGTCAATGAGGTGGACATGGAAGCAACGAACAGGAT
CAGGACCAGTGGATGATC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235782
Blue=ORF Red=Cloning site Green=Tag(s)

```
MSTLFPSLFPRVTETLWFNLDRPCVEETELQQEQQHQAQLQSI AEKDNLVPIGKPASEHYDDEEEED
DEDEDESEEDSEDEDEDMQDMDEMNDYNE SPDDGEVNEVDMEGNEQDQDQWMI
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRS NATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001278489

ORF Size: 363 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001278489.2](#)

RefSeq Size: 852 bp

RefSeq ORF: 366 bp

Locus ID: 25906


UniProt ID: [P60006](#)

Cytogenetics: 11q13.4

MW: 14.3 kDa

Gene Summary:

Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. In the complex, plays a role in the release of the mitotic checkpoint complex (MCC) from the APC/C: not required for APC/C activity itself, but promotes the turnover of CDC20 and MCC on the APC/C, thereby participating in the responsiveness of the spindle assembly checkpoint. Also required for degradation of CDC20.[UniProtKB/Swiss-Prot Function]

Product images:

Circular map for RG235782