

## Product datasheet for **RG235823**

### **C21orf7 (MAP3K7CL) (NM\_001286623) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C21orf7 (MAP3K7CL) (NM\_001286623) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** C21orf7  
**Synonyms:** C21orf7; HC21ORF7; TAK1L; TAKL; TAKL-1; TAKL-2; TAKL-4  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG235823 representing NM\_001286623.  
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTCAGCTGATTGCACCTTTAGAAGTTATGTGGAACGAGGCAGCAGATCTTAAGCCCCTTGCTCTG
TCACGCAGGCTGGAATGCAGTGGTGAATCATGGCTCACTACAGCCCTGACCTCTGGGCCAGAGATG
GAGTCTCGCTATTTGCCAGGTTGGTCTTGAACACCTGGCTTCAAGCAGTCTCTGCTTTGGCTTC
TTGAAGTCTGGATTACAGTATTTAGTCTTATGCTCTGCAACAAGTTGGCCATGTTGGAGGACAAT
CCAAAGGTCAGCAAGTTGGCTACTGGCGATTGGATGCTCACTCTGAAGCCAAAGTCTATTACTGTGCC
GTGGAAATCCCCAGCTCCCCTCTGGGT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

**Protein Sequence:** >Peptide sequence encoded by RG235823  
Blue=ORF Red=Cloning site Green=Tag(s)

```
MVQLIAPLEVMWNEADLKPLALSRRLECSGGIMAHYSPDLLPEMESRYFAQVGLLEHLASSPPAFGF
LKCLDYSISVLCSATSLAMLEDNPKVSKLATGDWMLTLKPKSITVPVEIPSSPLG
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPE
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001286623

**ORF Size:** 372 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**RefSeq:** [NM\\_001286623.2](#)

**RefSeq Size:** 1582 bp

**RefSeq ORF:** 375 bp


**Locus ID:** 56911

**Cytogenetics:** 21q21.3

**Protein Families:** Transcription Factors

**MW:** 13.8 kDa

Product images:



Circular map for RG235823