

Product datasheet for **RG235876**

PPCDC (NM_001301103) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PPCDC (NM_001301103) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PPCDC
Synonyms: coaC; MDS018; PPC-DC
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235876 representing NM_001301103.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAACCAAAGGCTCCTGTCCAGCTGCTGCACCCTTGATGGAGAGAAAATCCATGTTCTTGTGGGT
GTCACGGGGAGTGTGCGAGCCCTGAAGTTGCCTCTTCTGGTGTCAAAGCTTTTGGACATTCCTGGGACC
TGCGTCATGCGGGCCTGGGACCGCAGCAAGCCCTGCTCTTCTGCCGGCCATGAACACCGCCATGTGG
GAGCACCCGATCACAGCGCAGCAGGTAGACCAGCTCAAGGCCCTTGGCTATGTCGAGATCCCTGTGTG
GCCAAGAAGCTGGTGTGCGGAGATGAAGTCTCGGGCCATGGCTGAAGTGGGACCATCGTGGACAAA
GTGAAAGAAGTCTCTTCCAGCACAGTGGCTTCCAGCAGAGT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235876
Blue=ORF Red=Cloning site Green=Tag(s)

```
MEPKASCAPAAAPLMERKFHVLVGVTSVAALKLPLLVSKLLDIPGTCVMRAWDRSKPLLFCPAMNTAMW
EHPITAQQVDQLKAFGYVEIPCVAKKLVCGDEGLGAMAEVGTIVDKVKEVLFQHSQFQQS
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001301103

ORF Size: 387 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001301103.2](#)

RefSeq Size: 2043 bp

RefSeq ORF: 390 bp

Locus ID: 60490

UniProt ID: [Q96CD2](#)

Cytogenetics: 15q24.2


Protein Pathways: Metabolic pathways, Pantothenate and CoA biosynthesis

MW: 14.4 kDa

Gene Summary:

Biosynthesis of coenzyme A (CoA) from pantothenic acid (vitamin B5) is an essential universal pathway in prokaryotes and eukaryotes. PPCDC (EC 4.1.1.36), one of the last enzymes in this pathway, converts phosphopantothenoylcysteine to 4-prime-phosphopantetheine (Daugherty et al., 2002 [PubMed 11923312]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RG235876