

Product datasheet for **RG235895**

OPALIN (NM_001284322) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OPALIN (NM_001284322) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: OPALIN
Synonyms: HTMP10; TMEM10; TMP10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235895 representing NM_001284322.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGGATGACGTCTCTCTGTACAGGTGGGAAAGAAACGGACTGTGGCCCTCTCTGGATTAGCG
GCGGGCATAACCATTGCTGGTGGCCACAGCCCTGCTGGTGGCTTTACTATTTACTTTGATTCACCGAAGA
AGAAGCAGCATTGAGGCCATGGAGGAAAGTGACAGACCATGTGAAATTTAGAAATTGATGACAATCCC
AAGATATCTGAGAATCCTAGGAGATCACCCACACATGAGAAGAATACGATGGGAGCACAAGAGGCCAC
ATATATGTGAAGACTGTAGCAGGAAGCGAGGAACCTGTGCATGACCGTTACCGTCTACTATAGAAATG
GAAAGAAGGAGGGGATTGTGGTGGCTTGTGCCAGACTGAGCCTGGAA
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235895
Blue=ORF Red=Cloning site Green=Tag(s)

```
MRMTSSPVTGGKETDCGPSLGLAAGIPLLVATALLVALLFTLIHRRRSSIEAMEESDRPCEISEIDDNP
KISENPRRSPTHEKNTMGAQEAHIYVKTVAGSEEPVHRYRPTIEMERRRGLWWLVPRLSLE
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284322

ORF Size: 393 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001284322.1 , NP_001271251.1
RefSeq Size:	3717 bp
RefSeq ORF:	396 bp
Locus ID:	93377
UniProt ID:	Q96PE5
Cytogenetics:	10q24.1
Protein Families:	Transmembrane
MW:	15.1 kDa
Gene Summary:	Central nervous system-specific myelin protein that increase myelin genes expression during oligodendrocyte differentiation. Promotes oligodendrocyte terminal differentiation. [UniProtKB/Swiss-Prot Function]