

Product datasheet for **RG235902**

ZNF302 (NM_018675) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF302 (NM_018675) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZNF302
Synonyms: HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG235902 representing NM_018675.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCT
GCTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGGTCTTTCC
GTAACAAGCCATATGTGATCATGTTATTGGAGGATGAAAAGAGCCCTGGATGATGGAGAAAAAATCG
TCAAAAGCTTACCCATTTCTTTATCACACTCTGTTCTGCTTCTGTGAACCTTTGGATTCTCTGCTCTA
TTTGAGCATTGTCAGAAGTCACTGAAATATTTGAGTTGTCAGAACTATGTGTTTCTGGGTGCTTCAT
TTCTTATCCAATTCTCCTAATTCCACTGTAGAAGCTTTTTCAAGAAG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG235902
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MSQVTFSDVAIDFSHEEWACLDSAQRDL YKDVVMQNYENLVSVGLSVTKPYVIMLLEDGKEPMMMEKKL
SKAYPFPLSHSVPASVNFVGF SALFEHCSEVTEIFEL SELCVFWLHFLSNPNSTVEAFFKK
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: Sgfl-MluI



[View online »](#)


Cloning Scheme:


ACCN:	NM_018675
ORF Size:	393 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_018675.3
RefSeq Size:	3195 bp
RefSeq ORF:	396 bp
Locus ID:	55900
UniProt ID:	Q9NR11
Cytogenetics:	19q13.11
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	15.4 kDa

Gene Summary:

This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]

Product images:



Circular map for RG235902