

## Product datasheet for **RG236041**

### CCDC117 (NM\_001284265) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** CCDC117 (NM\_001284265) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** CCDC117  
**Synonyms:** dj366L4.1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG236041 representing NM\_001284265.  
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATCAGACAACGGAGAACCACAGTGTGAAGTTGCCGAAGGAAGCTTCAGGAGATTGAGGACAGG
ATAATTGATGAAGATGAAGAAGTTGAAGCTGACAGAAATGTTAACCATCTCCCGAGTCTTGTCTTTCT
GATACCATGAAAACAGGTTTGAAGAGGGAATTTGATGAAGTTTTTACAAAGAAAATGATTGAGTCTATG
AGCCGCTCTTCCATGGAGCTTGTCTCTGGAAACCCCTCCCTGAAGTCTTTCTGATAAGCCAAAGCCA
TCCTCTAATACTAAGAACTATACAGGAGAGAGCCAAGCTAAGCATGTAGCTGCTGGCACTGCCTTCCT
CAGAGAACTGAAGTGTTCGGAACCTCGGCCAACAGGGATGTCTTTTATAATAGTTTGGAGACAGCT
ACTAGCACAGAAGAAGAGATGGAAGTACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

**Protein Sequence:** >Peptide sequence encoded by RG236041  
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MDQTTGEPQCEVARRKLQEIEDRIIDEDEEVEADRNVNHLPSLVLSDTMKTGLKREFDEVFTKKMIESM
SRPSMELVLWKPLPELLSDKPKSSNTKNYTGESQAKHVAAGTAFPQRTELFSEPRPTGMSLYNSLETA
TSTEEEMEL
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_001284265

**ORF Size:** 441 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**RefSeq:** [NM\\_001284265.1](#), [NP\\_001271194.1](#)

**RefSeq Size:** 3717 bp

**RefSeq ORF:** 444 bp


**Locus ID:** 150275

**UniProt ID:** [Q8IWD4](#)

**Cytogenetics:** 22q12.1

**MW:** 17.2 kDa

Product images:



Circular map for RG236041