

Product datasheet for **RG236118**

Sigma1 receptor (SIGMAR1) (NM_001282205) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sigma1 receptor (SIGMAR1) (NM_001282205) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SIGMAR1
Synonyms: ALS16; DSMA2; hSigmaR1; OPRS1; SIG-1R; sigma1R; SR-BP; SR-BP1; SRBP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG236118 representing NM_001282205.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGTGGGCCGTGGCCGGCGGTGGCGTGGCCGCGCTGCTCCTGGCTGTCGCAGCGGTGCTGACC
CAGGTCGCTGGCTCTGGCTGGGTACGCAGAGCTTCGTCTCCAGCGCAAGAGATAGCGCAGTTGGCG
CGGCAGTACGCTGGCTGGACCACGAGCTGGCCTTCTCTCGTCTGATCGTGGAGCTGCGCGGCTGCAC
CCAGGCCACGTGCTGCCGACGAGGAGCTGCAGTGGGTGTTTCGTGAATGCGGGTGGCTGGATGGGCGCC
ATGTGCCTTCTGCACGCCTCGTGTCCGAGTATGTGCTGCTCTTCGGCACCGCCTGGGCTCCCGCGGC
CACTCGGGGCGCTACTGGGCTGAGATCTCGGATACCATCATCTCTGGCACCTTCCACCAGTGGAGAGAG
GGCACCACAAAAGTGAAGTCTTCTACCCAGGACCCCTTGACCAGCCAGGCC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG236118
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MQWAVGRRWAWAALLLAVA AVL TQVVWLWLTQSFVVFQREEIAQLARQYAGLDHELAFSRLIVELRRLH
PGHVL PDEELQWVFNAGGWMGAMCLLHASLSEYVLLFGTALGSRGHSGRYWAEISDTIISGTFHQWRE
GTTKSEVFYPGPLTSQA
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFSAIHP SILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282205

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001282205.2
RefSeq Size:	1510 bp
RefSeq ORF:	468 bp
Locus ID:	10280
UniProt ID:	Q99720
Cytogenetics:	9p13.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	17.9 kDa
Gene Summary:	This gene encodes a receptor protein that interacts with a variety of psychotomimetic drugs, including cocaine and amphetamines. The receptor is believed to play an important role in the cellular functions of various tissues associated with the endocrine, immune, and nervous systems. As indicated by its previous name, opioid receptor sigma 1 (OPRS1), the product of this gene was erroneously thought to function as an opioid receptor; it is now thought to be a non-opioid receptor. Mutations in this gene has been associated with juvenile amyotrophic lateral sclerosis 16. Alternative splicing of this gene results in transcript variants encoding distinct isoforms. [provided by RefSeq, Aug 2013]