

Product datasheet for **RG236293**

DCTN3 (NM_001281427) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DCTN3 (NM_001281427) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: DCTN3
Synonyms: DCTN-22; DCTN22
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG236293 representing NM_001281427.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGTCTGACTGACTTGCAGCGGTACAGGCCGAGTGGAAGAGCTGGAGCGCTGGGTGTACGGG
CCGGGCGGGCGCGCGGCTCACGGAAGGTGGTGACGGCCTGGTCAAGGTGCAGGTGGCTTTGGGGAAC
ATTTCCAGCAAGAGGGAGAGGGTGAAGATTCTCTACAAAAGATTGAAGATCTGATCAAGTACCTGGAT
CCTGAGTACATCGACCGCATTGCCATACCTGATGCCTCTAAGCTGCAATTCATCCTAGCAGAGGAGCAG
TTTATCCTTTCCAGTTGCACTCCTGGAGCAGGTGAATGCCTTGGTCCCATGCTGGACAGTGCTCAC
ATCAAAGCCGTTCTGAGCATGTGCCCCCTGCAGCGCTTGGCCAGATCCACATTCAGCAGCAGATT
CCTATTCTCACCTTGACCTACAACCAGACAATGCTTCTCTCCAAGCAATTCGTGCAGTGGGATGAGCTA
CTTTGCCAGCTAGAGGCCGCCACGCAAGTGAAGCCAGCAGAGGAG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

Protein Sequence: >Peptide sequence encoded by RG236293
Blue=ORF Red=Cloning site Green=Tag(s)

```
MAGLTDLQRLQARVEELERWVYGPGGARGSRKVADGLVKVQVALGNISSKRERVKILYKKIEDLIKYL
PEYIDRIAIPDASKLQFILAEEQFILSQVALLEQVNALVPMLDSAHIKAVPEHAARLQRLAQIHIQQQI
PILTLTYNQTMLLSKQFVQWDELLCQLEAATQVKPAEE
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

Restriction Sites: Sgfl-MluI



[View online >](#)

Cloning Scheme:


ACCN: NM_001281427

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001281427.2](#)

RefSeq Size: 827 bp

RefSeq ORF: 531 bp


Locus ID: 11258

Cytogenetics: 9p13.3

MW: 20.4 kDa

Gene Summary:

This gene encodes the smallest subunit of dynactin, a macromolecular complex consisting of 10 subunits ranging in size from 22 to 150 kD. Dynactin binds to both microtubules and cytoplasmic dynein. It is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle formation, cytokinesis, chromosome movement, nuclear positioning, and axonogenesis. This subunit, like most other dynactin subunits, exists only as a part of the dynactin complex. It is primarily an alpha-helical protein with very little coiled coil, and binds directly to the largest subunit (p150) of dynactin. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Product images:


Circular map for RG236293