

Cloning Scheme:


ACCN: NM_001288568

ORF Size: 549 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001288568.1](#), [NP_001275497.1](#)

RefSeq Size: 1743 bp

RefSeq ORF: 552 bp

Locus ID: 252839

UniProt ID: [Q9P0T7](#)


Cytogenetics: 1q32.1

Protein Families: Transmembrane

MW: 20.6 kDa

Gene Summary: May be involved in intracellular transport.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG236358