

## Product datasheet for **RG236400**

### ALS2CR15 (ICA1L) (NM\_001288624) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ALS2CR15 (ICA1L) (NM\_001288624) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ALS2CR15  
**Synonyms:** ALS2CR14; ALS2CR15  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG236400 representing NM\_001288624.  
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATTCCTTTGGGCAACCCAGACCAGAAGATAATCAGTCAGTAGTCAGAAGAATGCAAAAAGAAATAC
TGGAAAATAACAGGTCTTTATCAAAGCAACAGGAAAAAAGAGGATGAGCACTTGGTGGCGTCTGAT
GCTGAAGTGGATGCTAACTTGAGGTTTTTCACTCTGTTCAAGAGACATGCACTGAACCTTCTGAAGATA
ATCGAGAAATACCAGCTAAGACTCAATGTTATATCAGAGGAAGAAAATGAGCTAGGGCTCTTTTTAAAA
TTTCAAGCAGAACGGGATGCAACTCAAGCTGGCAAAATGATGGATGCCACTGGCAAGGCACCTTTGTCT
TCAGCCAAGCAAAGATTGGCCCTGTGTACTCCTCTGTCTCGTCTGAAGCAAGAAGTAGCAACATTCAGT
CAAAGGGCAGTATCTGATACCTTGATGACAATTAATCGGATGGAGCAGGCACGCACAGAATACAGAGGA
GCTCTACTGTGGATGAAAGATGTATCCCAAGAGCTGGACCCAGACACCTTAAAGCAAATGAAAAAGTTT
AGAAAAGTA
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

**Protein Sequence:** >Peptide sequence encoded by RG236400  
Blue=ORF Red=Cloning site Green=Tag(s)

```
MDSFGQRPEDNQSVVRRMQKKYWTKQVFIKATGKKEDEHLVASDAELDAKLEVFHSVQETCTELLKI
IEKYQLRLNVISEEENELGLFLKFAERDATQAGKMMDATGKALCSSAKQRLALCTPLSRLKQEVATFS
QRAVSDTLMTINRMEQARTEYRGALLWMKDVSQELDPDTLKQMEKFRKV
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPISILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```


**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


- ACCN:** NM\_001288624
- ORF Size:** 561 bp
- OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- RefSeq:** [NM\\_001288624.1](#), [NP\\_001275553.1](#)
- RefSeq Size:** 911 bp
- RefSeq ORF:** 564 bp
- Locus ID:** 130026
- UniProt ID:** [Q8NDH6](#)
- Cytogenetics:** 2q33.2
- MW:** 22 kDa

Product images:



Circular map for RG236400