

## Product datasheet for **RG236446**

### RHOH (NM\_001278368) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** RHOH (NM\_001278368) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** RHOH  
**Synonyms:** ARHH; TTF  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG236446 representing NM\_001278368.  
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGAGTTCCATCAAGTGCGTTGGTGGGCGACTCTGCTGTGGGAAAACCTCTCTGTTGGTGCGC
TTCACCTCCGAGACCTTCCCGGAGGCCTACAAGCCACAGTGTACGAGAACACAGGGGTGGACGCTTTC
ATGGATGGCATCCAGATCAGCCTGGCCTCTGGGACACAGCCGGCAATGACGCCTCAGAAGCATCCGG
CCCCTGTCCTACCAGCAGGCAGAGCTGGTGTGATGTGCTACTCTGTGGCAACCATAACTATTCTG
AACTTGAAGAACAAGTGGATTGGTGAAATTAGGAGCAACTTGCCCTGTACCCCTGTGCTGGTGGTGGCC
ACCCAGACTGACCAGCGGAGATGGGGCCCCACAGGGCCTCTGCGTCAATGCCATGGAAGGAAGAAA
CTGGCCCAGGATGTCAGAGCCAAGGGCTACCTGGAGTGCTCAGCCCTTAGCAATCGGGGAGTACAGCAG
GTGTTTGAGTGCCGCGTCCGAACTGCCTCAACCAGGCCAGGAGACGAAACAGAAGGAGGCTCTTCTCC
ATCAATGAGTGAAGATCTTC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```

**Protein Sequence:** >Peptide sequence encoded by RG236446  
Blue=ORF Red=Cloning site Green=Tag(s)

```
MLSSIKCVLVGDSAVGKTSLLVRFVTFPEAYKPTVYENTGVDVFMDCIQLSLGLWDTAGNDAFRSIR
PLSYQQADVLMCYSVANHNSFLNKNKWIIEIRSNLPTPVLVATQTDQREMGPFRASCVNAMEGKK
LAQDVRAGKYLECSALSNRGVQVFECAVRTAVNQARRRNRRLFSINECKIF
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```

**Restriction Sites:** Sgfl-MluI




**Cloning Scheme:**


- ACCN:** NM\_001278368
- ORF Size:** 573 bp
- OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- RefSeq:** [NM\\_001278368.2](#)
- RefSeq Size:** 2036 bp
- RefSeq ORF:** 576 bp
- Locus ID:** 399
- UniProt ID:** [Q15669](#)
- Cytogenetics:** 4p14
- Protein Families:** Transcription Factors
- Protein Pathways:** Leukocyte transendothelial migration
- MW:** 21.3 kDa

**Gene Summary:**

The protein encoded by this gene is a member of the Ras superfamily of guanosine triphosphate (GTP)-metabolizing enzymes. The encoded protein is expressed in hematopoietic cells, where it functions as a negative regulator of cell growth and survival. This gene may be hypermutated or misexpressed in leukemias and lymphomas. Chromosomal translocations in non-Hodgkin's lymphoma occur between this locus and B-cell CLL/lymphoma 6 (BCL6) on chromosome 3, leading to the production of fusion transcripts. Alternative splicing in the 5' untranslated region results in multiple transcript variants that encode the same protein. [provided by RefSeq, May 2013]

**Product images:**

Circular map for RG236446