

Product datasheet for **RG236696**

METTL11A (NTMT1) (NM_001286797) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	METTL11A (NTMT1) (NM_001286797) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	METTL11A
Synonyms:	AD-003; C9orf32; HOMT1A; METTL11A; NRMT; NRMT1; NTM1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG236696 representing NM_001286797. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACGAGCGAGGTGATAGAAGACGAGAAGCAATTCTATTCCAAGGCCAAGACCTACTGGAAACAAATC
CCACCCACGGTGGACGGCATGCTTGGGGGGTATGGCCACATCTCCAGCATCGACATCAACAGCTCCCGG
AAGTTTCTGCAGAGGTTTTTGGGGAAGGCCCAACAAGACAGGAACGTCTGTGCCCTGGACTGTGGA
GCTGGCATTGGGAGGATACCAAGCGGCTGCTCCTGCCGCTGTTCCAGAGAGGTGGATATGGTCGACATA
ACGGAGGACTTCCTGGTCAAGCCAAGACCTACCTGGGGGAGGAGGCAAGAGGGTGAGGAACTACTTC
TGTTGTGGGCTCCAGGACTTCACCCCGGAGCCGACTCTTACGACGTGATCTGGATCCAGTGGGTGATA
GGCCACCTCACCGATCAGCACCTGGCCGAGTTCCTGCGGCGCTGCAAGGGCAGCCTCCGCCCAACGGC
ATCATCGTCATCAAAGACAACATGGCCCAGGAGGGCGTGATTCTGGACGACGTGGACAGCAGCGTGTGC
CGGGACCTTGACGTGGTCCGCAGGATCATCTGCAGTGCAGGCCTCAGCCTCCTGGCCGAGGAGGGCAG
GAGAACCTCCCGATGAGATCTACCATGTCTATAGCTTTGCCCTGAGA
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG236696
 Blue=ORF Red=Cloning site Green=Tag(s)

MTSEVIEDEKQFYSAKTYWKQIPPTVDGMLGGYGHISSIDINSSRKFLQRFLREGPNKGTSCALDCG
 AGIGRITKRLLLPLFREVDMDITEDFLVQAKTYLGEEGRVRYFCCLQDFTPEPDSYDVIWIQWVI
 GHLDQHLAEFLRRCKGSLRPNGIIVIKDNMAQEGVILDDVDSSVCRDLDVVRIICSAGLSLLAEERQ
 ENLPDEIYHVVSFALR
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIIGDFKVMGTGFPEP
 SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001286797

ORF Size: 669 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


RefSeq: [NM_001286797.1](#), [NP_001273726.1](#)

RefSeq Size: 1501 bp

RefSeq ORF: 672 bp

Locus ID:	28989
UniProt ID:	Q9BV86
Cytogenetics:	9q34.11
Protein Families:	Druggable Genome
MW:	25.4 kDa
Gene Summary:	The METTL11A gene encodes an N-terminal methyltransferase for the RAN (MIM 601179) guanine nucleotide exchange factor regulator of chromosome condensation 1 (RCC1; MIM 179710). METTL11A enzyme alpha-N-methylates other protein targets such as SET (MIM 600960) and RB (MIM 180200).[supplied by OMIM, Nov 2010]

Product images:



Circular map for RG236696