

Product datasheet for **RG236813**

CD40 (NM_001302753) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD40 (NM_001302753) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CD40
Synonyms:	Bp50; CDW40; p50; TNFRSF5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG236813 representing NM_001302753. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTCGTCTGCCTCTGCAGTGCCTCTGGGGCTGCTTGCTGACCGCTGTCCATCCAGAACCACCC
ACTGCATGCAGAGAAAAACAGTACCTAATAAACAGTCAGTGCTGTTCTTTGTGCCAGCCAGGACAGAAA
CTGGTGTGACTGCACAGAGTTCAGTAAACGGAATGCCTTCCTTGGGTGAAAGCGAATTCCTAGAC
ACCTGGAACAGAGAGACACTGCCACCAGCACAATACTGCGACCCCAACCTAGGGCTTCGGGTCCAG
CAGAAGGGCACCTCAGAAACAGACACCATCTGCACCTGTGAAGAAGGCTGGCACTGTACGAGTGAGGCC
TGTGAGAGCTGTGCTCCTGCACCGCTCATGCTCGCCGGCTTTGGGGTCAAGCAGATTGCTACAGGGGTT
TCTGATACCATCTGCGAGCCCTGCCAGTCGGCTTCTTCTCCAATGTGTCATCTGCTTTCGAAAAATGT
CACCCCTTGGACAAGCTGTGAGACCAAAGACCTGTTGTGCAACAGGCAGGCACAAACAAGACTGATGTT
GTCTGTGGTGAGTCTGGACAATGGGCCCTGGAGAAAGCCTAGGAAGGTCCCGGATCGGCTGAGAGC
CCTGGTGGTGATCCCATCATCTTCGGGATCCTGTTTGCCATCCTCTTGGTGCTGGTCTTTATCAAAAA
GGTGGCCAAGAAGCCAACCAA
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG236813
 Blue=ORF Red=Cloning site Green=Tag(s)


MVRLPLQCVLWGCLLTAVHPEPPTACREKQYLINSQCSSLCPGQKLVSDCTEFTETECLPCGESEFLD
 TWNRETHCHQHKYCDPNLGLRVQKGTSETDTICTCEEGWHCTSEACESCVLHRS CSPGFGVKQIATGV
 SDTICEPCVGVFFSNVSSAFEKCHPWTSCETKDLVVQAGTNKTDVVCGESWTMGPGESLGRSPGSAES
 PGGDPHHLRDPVCHPLGAGLYQKGGQEANQ
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSGYENPFLHAINNGGYTNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
 SVIFTDKIIRS NATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001302753

ORF Size:	711 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001302753.2
RefSeq Size:	1669 bp
RefSeq ORF:	714 bp
Locus ID:	958
UniProt ID:	P25942
Cytogenetics:	20q13.12
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Allograft rejection, Asthma, Autoimmune thyroid disease, Cell adhesion molecules (CAMs), Cytokine-cytokine receptor interaction, Primary immunodeficiency, Systemic lupus erythematosus, Toll-like receptor signaling pathway, Viral myocarditis
MW:	26.2 kDa
Gene Summary:	This gene is a member of the TNF-receptor superfamily. The encoded protein is a receptor on antigen-presenting cells of the immune system and is essential for mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development, and germinal center formation. AT-hook transcription factor AKNA is reported to coordinately regulate the expression of this receptor and its ligand, which may be important for homotypic cell interactions. Adaptor protein TNFR2 interacts with this receptor and serves as a mediator of the signal transduction. The interaction of this receptor and its ligand is found to be necessary for amyloid-beta-induced microglial activation, and thus is thought to be an early event in Alzheimer disease pathogenesis. Mutations affecting this gene are the cause of autosomal recessive hyper-IgM immunodeficiency type 3 (HIGM3). Multiple alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported. [provided by RefSeq, Nov 2014]