

Product datasheet for **RG236861**

Chimaerin 2 (CHN2) (NM_001293075) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chimaerin 2 (CHN2) (NM_001293075) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CHN2
Synonyms:	ARHGAP3; BCH; CHN2-3; RHOGAP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG236861 representing NM_001293075. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTCTCTGAAGAACTGTGGCTGGAAAATGAGAAAAAGTGTGCTGTGGTTCGGAAGTCTAAGCAGGGC
AGGAAACGCCAAGAACTGCTGGCCGTAGCCTTCGGGGTGAAGGTCCACACGTTCCGAGGCCACACTGG
TGTGAATATTGTCCAATTTTCATGTGGGGGCTCATCGCCCAAGGGGTCCGGTCTCAGACTGTGGATTG
AACGTACACAAACAGTGTCCAAGCACGTTCCAATGACTGCCAACCTGATCTCAAGAGGATCAAGAAA
GTGTACTGTTGTGACCTCACAACACTTGTGAAGGCTCACAACACTCAGAGACCCATGGTGGTAGACATA
TGCATTCGGGAAATTGAAGCAAGAGGATTAATAATCGGAAGGCCTTACAGAGTCTCTGGGTTCACTGAA
CACATTGAAGATGTCAAATGGCATTGACAGAGAAATCTCCAATGCAGATGAGAGGCTGGAAGCCGTC
CATGAAGTGTGATGCTGCTGCCTCCTGCCACTATGAAACCCTCCGGTACCTAATGATCCACCTCAA
AAGGTTACTATGAATGAAAAAGACAATTTTCATGAATGCAGAAAATCTGGGGATCGTGTGGGCCCACT
CTGATGAGGCCCCCTGAGGACAGCACCTGACCACCTGCATGATATGCGGTACCAAAAGCTGATTGTG
CAGATTTAATAGAAAACGAAGACGTTTTATTC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG236861
 Blue=ORF Red=Cloning site Green=Tag(s)


MFSEELWLENEKKCAVVRKSKQGRKRQELLAVAFGVKVTFRGPHWCEYCANFMWGLIAQGVRCSDCGL
 NVHKQCSKHVPNDQCQDLKRIKKVYCCDLTTLVKAHNTQRPMVVDICIREIEARGLKSEGLYRVSGFTE
 HIEDVKMAFDREISNADERLEAVHEVLMMLPPAHYETLRYLMIHLKKVTMNEKDNFMNAENLGIVFGPT
 LMRPPEDSTLTLHDMRYQKLIQVILINEDVLF
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYTNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001293075

ORF Size:	723 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001293075.1 , NP_001280004.1
RefSeq Size:	2742 bp
RefSeq ORF:	726 bp
Locus ID:	1124
UniProt ID:	P52757
Cytogenetics:	7p14.3
MW:	28.4 kDa
Gene Summary:	This gene encodes a guanosine triphosphate (GTP)-metabolizing protein that contains a phorbol-ester/diacylglycerol (DAG)-type zinc finger, a Rho-GAP domain, and an SH2 domain. The encoded protein translocates from the cytosol to the Golgi apparatus membrane upon binding by diacylglycerol (DAG). Activity of this protein is important in cell proliferation and migration, and expression changes in this gene have been detected in cancers. A mutation in this gene has also been associated with schizophrenia in men. Alternative transcript splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, May 2014]