

Product datasheet for **RG237030**

AMD1 (NM_001287215) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AMD1 (NM_001287215) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AMD1
Synonyms:	ADOMETDC; AMD; SAMDC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237030 representing NM_001287215. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTGTCTCCAAGAGACGTTTCATTTTGAAGACATGTGGTACCACCCTCTTGCTGAAAGCACTGGTT
CCCCTGTTGAAGCTTGCTAGGGATTACAGTGGTTTGACTCAATCAAAGCTTCTTTATTCTCGTAAG
AATTTTCATGAAGCCTTCTACCAAGGGTACCCACACCGGAATTTCCAGGAAGAAATAGAGTTTCTTAAT
GCAATTTCCCAAATGGAGCAGCATATTGTATGGGACGTATGAATTCTGACTGTTGGTACTTATACT
CTGGATTTCCAGAGAGTCGGTAAATCAGTCAGCCAGATCAAACCTTGGAAATTCGATGAGTGAGCTT
GACCCAGCAGTTATGGACCAGTTCTACATGAAAGATGGTGTACTGCAAAGGATGTCACCTCGTGAGAGT
GGAATTCGTGACCTGATACCAGTTCTGTCAATGATGCCACAATGTTCAATCCTTGTGGGTATTCGATG
AATGGAATGAAATCGGATGGAACCTATTGGACTATTCACATCACTCCAGAACCAGAAATTTCTTATGTT
AGCTTTGAAACAACTTAAGTCAGACCTCCTATGATGACCTGATCAGGAAAGTTGTAGAAGTCTTCAAG
CCAGGAAAATTTGTGACCACCTTGTGTTGTAATCAGAGTTCTAAATGTCGCACAGTGCTTGTTCGCC
CAGAAGATTGAAGTTTTAAGCGTCTTGATTGCCAGAGTGCTATGTTCAATGATTACAATTTGTTTTT
ACCAGTTTTGCTAAGAAGCAGCAACAACAGCAGAGT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAAAC
```



[View online »](#)

ORF Size:	795 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001287215.1 , NP_001274144.1
RefSeq Size:	3200 bp
RefSeq ORF:	798 bp
Locus ID:	262
Cytogenetics:	6q21
Protein Families:	Druggable Genome
Protein Pathways:	Arginine and proline metabolism, Cysteine and methionine metabolism, Metabolic pathways
MW:	30.9 kDa
Gene Summary:	This gene encodes an important intermediate enzyme in polyamine biosynthesis. The polyamines spermine, spermidine, and putrescine are low-molecular-weight aliphatic amines essential for cellular proliferation and tumor promotion. Multiple alternatively spliced transcript variants have been identified. Pseudogenes of this gene are found on chromosomes 5, 6, 10, X and Y. [provided by RefSeq, Dec 2013]