

Product datasheet for **RG237138**

PPP1R7 (NM_001282410) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP1R7 (NM_001282410) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PPP1R7
Synonyms:	SDS22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237138 representing NM_001282410. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCGGAACGCGCGCGGGGCAGCAACAGTCGCAGGAGATGATGGAGGTTGACAGGCGGGTCGAG
TCTGAAGAATCCGGCGATGAAGAAGGGAAGAAACACAGCAGTGGCATCGTGGCCGACCTCAGTGAACAG
AGCCTGAAGGATGGGAGGAGCGGGGGAGGAGGACCCAGAAGAAGAATGAGCTGCCTGTGGACATG
GAAACCATCAACCTGGACAGAGATGCAGAGGATGTTGATTTGAATCACTATCGCATAGGGAAGATTGAA
GGATTTGAGGTACTGAAGAAAGTGAAGACTCTCTGCCTCCGCCAAAATTTAATTAATGCATTGAGAAT
CTGGAGGAGCTACAGAGTCTTCGAGAGCTGGATCTTTACGACAACCAGATCAAGAAGATTGAGAATCTG
GAGGCGCTAACAGAGCTGGAGATTCTAGATATTTCTTTAATCTGCTGAGAAACATCGAAGGGGTTGAC
AAGTTGACACGACTGAAAAACTCTTCTTGGTCAACAATAAAATCAGTAAAATTGAGAACTTAAGCAAC
TTACATCAACTACAGATGCTAGAGCTGGGATCTAACCGCATCCGGGCAATCGAAAATATCGACACCTTA
ACCAACCTGGAGAGTTTGTTTTTGGGAAAAACAAAATTAATAAATTCAGAACCTGGATGCGCTCACC
AACCTGACAGTCTCAGTATGCAGAGCAACCGGCTGACCAAGATCGAGGGTCTGCAGAACCTGGTGAAAC
CTGCGGGAGCTGTACCTTAGCCACAATGGCATCGAGGTCATCGAGGGCCTGGAGAACAATGTGCAGGAC
AGCCTCACGTAC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237138
 Blue=ORF Red=Cloning site Green=Tag(s)

MAAERGAGQQSQEMMEVDRRVESEESGDEEGKKHSSGI VADLSEQSLKDGEERGEEDPEEEHEL PVDM
 ETINLDRDAEDVDLNHYRIGKIEGFVLLKVKTLCLRQNL IKCIENLEELQSLRELDLYDNQIKKIENL
 EALTELEILDISFNLLRNIEGVDKLTRLKKLFLVNNKISKIENLSNLHQLQMLELGSNRIRAIENIDTL
 TNLESFLGKNKITKLQNL DALTNLTVLSMQSNRLTKIEGLQNLVNLRELYL SHNGIEVIEGLENNVQD
 SLTY
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRS NATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001282410

ORF Size: 840 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


RefSeq: [NM_001282410.1](#), [NP_001269339.1](#)

RefSeq Size: 1803 bp

RefSeq ORF: 843 bp

Locus ID: 5510
UniProt ID: [Q15435](#)
Cytogenetics: 2q37.3
Protein Families: Druggable Genome, Phosphatase
MW: 32.5 kDa
Gene Summary: This gene encodes a protein subunit that regulates the activity of the serine/threonine phosphatase, protein phosphatase-1. The encoded protein is required for completion of the mitotic cycle and for targeting protein phosphatase-1 to mitotic kinetochores. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]

Product images:



Circular map for RG237138