

Product datasheet for **RG237141**

LAIR1 (NM_001289026) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LAIR1 (NM_001289026) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LAIR1
Synonyms:	CD305; LAIR-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237141 representing NM_001289026. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAAGAGAAATGGTCTCTGCCTGGCCAGACCATCCACACGCAGGAGGATCTGCCAGACCTCC
ATCTCGGCTGAGCCAGGCACCGTGATCCCCCTGGGAGCCATGTGACTTTCGTGTGCCGGGGCCCGTT
GGGGTTCAAACATTCGCCTGGAGAGGGACAGTAGATCCACATACAATGATACTGAAGATGTGTCTCAA
GCTAGTCCATCTGAGTCAGAGGCCAGATCCGCATTGACTCAGTAAGAGAAGGAAATGCCGGCTTTAT
CGCTGCATCTATTATAAGCCCCCTAAATGGTCTGAGCAGAGTGACTACCTGGAGCTGCTGGTGAAGAA
AGCTCTGGAGGCCCGGACTCCCCGGACACAGAGCCGGCTCCTCAGCTGGACCCACGCAGAGGCCGTCG
GACAACAGTCACAATGAGCATGCACCTGCTTCCAAGGCCTGAAAGCTGAGCATCTGTATATTCTCATC
GGGTCTCAGTGGTCTTCTCTTCTGTCTCCTCCTCCTGGTCTTCTGCTCCATCGCCAGAATCAG
ATAAAGCAGGGGCCCCCAGAAGCAAGGACGAGGAGCAGAAGCCACAGCAGAGGCCTGACCTGGCTGTT
GATGTTCTAGAGAGGACAGCAGACAAGGCCACAGTCAATGGACTTCTGAGAAGGACAGAGAGACGGAC
ACCTCGGCCCTGGCTGCAGGGAGTTCCAGGAGGTGACGTATGCTCAGCTGGACCACTGGGCCCTCACA
CAGAGGACAGCCCGGCTGTGTCCCACAGTCCACAAAGCCATGGCCGAGTCCATCACGTATGCAGCC
GTTGCCAGACAC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237141
 Blue=ORF Red=Cloning site Green=Tag(s)

MEREMVLCLAQTIHTQEDLPRPSISAEPGTVIPLGSHVTFVCRGPVGVQTFRLERDSRSTYNDTEDVSQ
 ASPSESEARFRIDSVREGNAGLYRCIYYKPPKWEQSDYLELLVKESGGPDSPDEPGSSAGPTQRPS
 DNSHNEHAPASQGLKAEHLYILIGVSVVFLFCLLLLVLFLHRQNIKQGPFRSKDEEQKPQQRDLAV
 DVLERTADKATVNLPEKDRETDTSALAAGSSQEVTYAQLDHWALTQRTARAVSPQSTKPMAESITYAA
 VARH
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
 SVIFTDKIIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001289026

ORF Size: 840 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001289026.2](#), [NP_001275955.2](#)


RefSeq Size: 2766 bp

RefSeq ORF: 843 bp

Locus ID: 3903
Cytogenetics: 19q13.42
Protein Families: Transmembrane
MW: 31.3 kDa

Gene Summary: The protein encoded by this gene is an inhibitory receptor found on peripheral mononuclear cells, including natural killer cells, T cells, and B cells. Inhibitory receptors regulate the immune response to prevent lysis of cells recognized as self. The gene is a member of both the immunoglobulin superfamily and the leukocyte-associated inhibitory receptor family. The gene maps to a region of 19q13.4 called the leukocyte receptor cluster, which contains at least 29 genes encoding leukocyte-expressed receptors of the immunoglobulin superfamily. The encoded protein has been identified as an anchor for tyrosine phosphatase SHP-1, and may induce cell death in myeloid leukemias. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Product images:



Circular map for RG237141