

Product datasheet for **RG237201**

LAIR1 (NM_001289025) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LAIR1 (NM_001289025) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LAIR1
Synonyms:	CD305; LAIR-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237201 representing NM_001289025. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTCCCCACCCACCGCCCTCTGGGCTAGTGCTCTGCCTGGCCAGACCATCCACACGCAGGAG
GATCTGCCAGACCCCTCCATCTCGGCTGAGCCAGGCACCGTGATCCCCCTGGGGAGCCATGTGACTTTC
GTGTGCCGGGGCCGGTTGGGGTTCAAACATTCGCCTGGAGAGGGACAGTAGATCCACATAAATGAT
ACTGAAGATGTGTCTCAAGCTAGTCCATCTGAGTCAGAGGCCAGATTCCGCATTGACTCAGTAAGAGAA
GGAAATGCCGGGCTTTATCGCTGCATCTATTATAAGCCCCCTAAATGGTCTGAGCAGAGTGACTACCTG
GAGCTGTGGTAAAAGAAAGCTCTGGAGCCCGGACTCCCGGACACAGAGCCCGCTCCTCAGCTGGA
CCCACGCAGAGGCCGTGCGACAACAGTCACAATGAGCATGCACCTGCTTCCCAAGGCTGAAAGCTGAG
CATCTGTATATTCTCATCGGGTCTCAGTGGTCTTCTCTTCTGTCTCCTCCTCCTGTCCTCTTCTGC
CTCCATCGCCAGAATCAGATAAAGCAGGGGCCCCCAGAAGCAAGGACGAGGAGCAGAAGCCACAGCAG
AGGCTGACCTGGCTGTTGATGTTCTAGAGAGGACAGCAGACAAGGCCACAGTCAATGGACTTCTGAG
AAGGACAGAGAGACGGACACCTCGGCCCTGGCTGCAGGGAGTCCAGGAGGTGACGTATGCTCAGCTG
GACCACTGGGCCCTCACACAGAGGACAGCCCGGGCTGTGTCCCCACAGTCCACAAAGCCCATGGCCGAG
TCCATCACGTATGCAGCCGTTGCCAGACAC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237201
 Blue=ORF Red=Cloning site Green=Tag(s)

MSPHPTALLGLVLCLAQTIHTQEDLPRPSISAEPGTVIPLGSHVTFVCRGPVGVQTFRLERDSRSTYND
 TEDVVSQASPSESEARFRIDSVREGNAGLYRCIYKPPKWSEQSDYLELLVKESSGGPDSPDTEPGSSAG
 PTQRPSDNSHNEHAPASQGLKAEHLIYLIGVSVVFLFCLLLLVLFLCHRQNKIKQGPGRSKDEEQKPKQ
 RPDLAVDLERTADKATVNGLPEKDRETDTSALAAAGSSQEVTYAQLDHWALTQRTARAVSPQSTKPMAE
 SITYAAVARH
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001289025

ORF Size: 858 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001289025.3](#)


RefSeq Size: 2815 bp

RefSeq ORF: 861 bp

Locus ID: 3903
UniProt ID: [Q6GTX8](#)
Cytogenetics: 19q13.42
Protein Families: Transmembrane
MW: 31.8 kDa

Gene Summary: The protein encoded by this gene is an inhibitory receptor found on peripheral mononuclear cells, including natural killer cells, T cells, and B cells. Inhibitory receptors regulate the immune response to prevent lysis of cells recognized as self. The gene is a member of both the immunoglobulin superfamily and the leukocyte-associated inhibitory receptor family. The gene maps to a region of 19q13.4 called the leukocyte receptor cluster, which contains at least 29 genes encoding leukocyte-expressed receptors of the immunoglobulin superfamily. The encoded protein has been identified as an anchor for tyrosine phosphatase SHP-1, and may induce cell death in myeloid leukemias. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Product images:



Circular map for RG237201