

## Product datasheet for **RG237303**

### Mu Opioid Receptor (OPRM1) (NM\_001285526) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Mu Opioid Receptor (OPRM1) (NM_001285526) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mu Opioid Receptor
Synonyms:	LMOR; M-OR-1; MOP; MOR; MOR1; OPRM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237303 representing NM_001285526. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGAAGACTGCCACCAACATCTACATTTTCAACCTTGCTCTGGCAGATGCCTTAGCCACCAGTACCCCTG
CCCTTCCAGAGTGTGAATTACCTAATGGGAACATGGCCATTTGGAACCATCCTTTGCAAGATAGTGATC
TCCATAGATTACTATAACATGTTCAACAGCATATTCACCCTCTGCACCATGAGTGTGATCGATACATT
GCAGTCTGCCACCCTGTCAAGGCCTTAGATTTCCGTA CTCCCGAAATGCCAAAATTATCAATGTCTGC
AACTGGATCCTCTCTTCAGCCATTGGTCTTCTGT AATGTT CATGGCTACAACAAAATACAGGCAAGGT
TCCATAGATTGTACACTAACATTCTCTCATCCAACCTGGTACTGGGAAAACCTGCTGAAGATCTGTGTT
TTCATCTTCGCCTTCATTATGCCAGTGCTCATATTACCGTGTGCTATGGACTGATGATCTTGGCCCTC
AAGAGTGTCCGCATGCTCTCTGGCTCCAAGAAAAGGACAGGAATCTTCGAAGGATCACCAGGATGGTG
CTGGTGGTGGTGGCTGTGTTTCATCGTCTGCTGGACTCCCATTACATTTACGTCATCATTAAAGCCTTG
GTTACAATCCCAGAAACTACGTTCCAGACTGTTTCTTGGCACTTCTGCATTGCTCTAGGTTACACAAAC
AGCTGCCTCAACCCAGTCCTTTATGCATTTCTGGATGAAAACCTCAAACGATGCTTCAGAGAGTTCTGT
ATCCCAACCTCTCCAACATTGAGCAACAAAACCTCACTCGAATTCGTCAGAACACTAGAGACCACCCC
TCCACGGCCAATACAGTGGATAGAATAATCATCAGCTAGAAAATCTGGAAGCAGAAACTGCTCCGTTG
CCC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

**Protein Sequence:** >Peptide sequence encoded by RG237303  
 Blue=ORF Red=Cloning site Green=Tag(s)

MKTATNIYIFNLALADALATSTLPPFQSVNYLMGTWPFGTILCKIVISIDYNNMFTSIFTLCTMSVDRII  
 AVCHPVKALDFRTPRNAKIINVCNWILSSAIGLPVMFMATTKYRQGSIDCTLTFSHPTWYENLLKICV  
 FIFAFIMPVLIITVCYGLMILRLKSVRMLSGSKEKDRNLRRITRMVLVVAVFIVCWTPIHIVYIKAL  
 VTIPETTFQTVSWHFCIALGYTNSCLNPVLYAFLDENFKRCFREFCRIPTSSNIEQQNSTRIRQNRDHP  
 STANTVDRTNHQLLENLEAETAPLP  
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV  
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED  
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA  
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001285526

**ORF Size:** 900 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**RefSeq:** [NM\\_001285526.1](#), [NP\\_001272455.1](#)

**RefSeq Size:** 15825 bp


**RefSeq ORF:** 903 bp

**Locus ID:** 4988  
**UniProt ID:** [P35372](#)  
**Cytogenetics:** 6q25.2  
**Protein Families:** Druggable Genome, GPCR, Transmembrane  
**Protein Pathways:** Neuroactive ligand-receptor interaction

**MW:** 34.4 kDa

**Gene Summary:** This gene encodes one of at least three opioid receptors in humans; the mu opioid receptor (MOR). The MOR is the principal target of endogenous opioid peptides and opioid analgesic agents such as beta-endorphin and enkephalins. The MOR also has an important role in dependence to other drugs of abuse, such as nicotine, cocaine, and alcohol via its modulation of the dopamine system. The NM\_001008503.2:c.118A>G allele has been associated with opioid and alcohol addiction and variations in pain sensitivity but evidence for it having a causal role is conflicting. Multiple transcript variants encoding different isoforms have been found for this gene. Though the canonical MOR belongs to the superfamily of 7-transmembrane-spanning G-protein-coupled receptors some isoforms of this gene have only 6 transmembrane domains. [provided by RefSeq, Oct 2013]

### Product images:



Circular map for RG237303