

Product datasheet for **RG237331**

WDR61 (NM_001303247) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR61 (NM_001303247) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	WDR61
Synonyms:	REC14; SKI8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237331 representing NM_001303247. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGACCAACCAGTACGGTATTCTCTTCAAACAAGAGCAAGCCCATGATGATGCCATTTGGTCAGTTGCT
TGGGGGACAAACAAGAAGGAAAACCTCTGAGACAGTGGTCACAGGCTCCCTAGATGACCTGGTGAAGGTC
TGGAAATGGCGTGATGAGAGGCTGGACCTACAGTGGAGTCTGGAGGGACATCAGCTGGGAGTGGTGTCT
GTGGACATCAGCCACACCCTGCCATTGCTGCATCCAGCTCTCTTGATGCTCATATTCGCTTTGGGAC
TTGGAAAATGGCAAACAGATAAAGTCCATAGATGCAGGACCTGTGGATGCCTGGACTTTGGCCTTTTCT
CCTGATCCCAGTATCTGGCCACAGGAACCTCATGTCGGGAAAGTGAACATTTTTGGTGTGAAAGTGGG
AAAAAGGAATATTCTTTGGACACGAGAGGAAAATTCATTCTTAGTATTGCATATAGTCCTGATGGGAAA
TACCTAGCCAGTGGAGCCATAGATGGAATCATCAATTTTTGATATTGCAACTGGAAAACCTTTCGCAT
ACCCTGGAAGGCCATGCCATGCCATTGCTCCTTGACCTTTTCCCGGACTCCCAGCTCCTTGTCACT
GCTTCAGATGATGGCTACATCAAGATCTATGATGTACAACATGCCAATTTGGCTGGCAGCGTGAGCGGC
CATGCCTCCTGGGTGCTGAACGTTGCATTCTGTCTGATGACACTCACTTTGTTCCAGTTGCTGTGAC
AAAAGTGTAAGTTTGGGATGTTGGAACGAGGACTTGTGTTACACCTCTTTGATCACCAGGATCAG
GTCTGGGAGTAAAATACAATGGAAATGGTTCAAAAATGTGTCTGTTGGAGATGACCAGGAAATTCAC
ATCTATGATTGTCCAATT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237331
 Blue=ORF Red=Cloning site Green=Tag(s)


MTNQYGILFKQEQAHDDAIWSVAWGNTNKKENSETVVTGSLDDLKVKVWWRDERLDLQWSLEGHQLGVVS
 VDISHTLPIAASSSLDAHIRLWLENGKQIKSIDAGPVDATLAFSPDSQYLATGTHVGKVNIFGVESG
 KKEYSLDRGKFI LSIAYSPDGKYLASGAIDGI INIFDIATGKLLHTLEGHAMP I RSLTFSPDSQLLV T
 ASDDGYIKIYDVQHANLAGT LSGHASWVLNVAFCPDDTHFVSSSSDKSVK VWDVGT R TCVHTFFDHQDQ
 VWGVKYNNGNSKIVSVGDDQEIHIYDCPI
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYTNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
 SVIFTDKII RSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001303247

ORF Size:	915 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001303247.2
RefSeq Size:	1439 bp
RefSeq ORF:	918 bp
Locus ID:	80349
UniProt ID:	Q9GZS3
Cytogenetics:	15q25.1
Protein Pathways:	RNA degradation
MW:	34 kDa
Gene Summary:	WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).[supplied by OMIM, Mar 2008]