

Product datasheet for **RG237524**

MSF (SEPT9) (NM_001293697) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MSF (SEPT9) (NM_001293697) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SEPTIN9
Synonyms:	AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237524 representing NM_001293697. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGCCGACACCCCCAGAGATGCCGGGCTCAAGCAGGCGCCTGCATCACGGAACGAGAAGGCCCCGGTG
GACTTCGGCTACGTGGGATTGACTCCATCCTGGAGCAGATGCGCCGGAAGGCCATGAAGCAGGGCTTC
GAGTTCAACATCATGGTGGTCGGGCAGAGCGGCTTGGTAAATCCACCTAATCAACACCCTCTTCAAA
TCCAAAATCAGCCGGAAGTCGGTGCAGCCACCTCAGAGGAGCGCATCCCAAGACCATCGAGATCAAG
TCCATCACGCACGATATTGAGGAGAAAGGCGTCCGGATGAAGCTGACAGTGATTGACACACCAGGGTTC
GGGACCACATCAACAACGAGAAGTCTGGCAGCCCATCATGAAGTTCATCAATGACCAGTACGAGAAA
TACCTGCAGGAGGAGGTCAACATCAACCGCAAGAAGCGCATCCCGGACACCCGCTCCACTGCTGCCTC
TACTTCATCCCCGCCACCGGCCACTCCCTCAGGCCCTGGACATCGAGTTTATGAAACGCCTGAGCAAG
GTGGTCAACATCGTCCCTGTCAATCGCAAGGCGGACACTCACCTGGAGGAGAGGGTCCACTTCAAA
CAGCGGATCACCGCAGACCTGCTGTCCAACGGCATCGACGTGTACCCCAAGGAATTTGATGAGGAC
TCGGAGGACCGGCTGGTGAACGAGAAGTCCGGGAGATGATCCCATTTGCTGTGGTGGCAGTGACCAC
GAGTACCAGGTCAACGGCAAGAGGATCCTTGGGAGGAAGCAAGTGGGTACCATCGAAGTTGAAAAC
ACCACACTGTGAGTTTGCCTACCTGCGGGACCTTTCATCAGGACGCACATGCAGAACATCAAGGAC
ATCACCAGCAGCATCCACTTCGAGGCGTACCGTGTGAAGCGCTCAACGAGGGCAGCAGCGCCATGGCC
AACGGCATGGAGGAGAAGGAGCCAGAAGCCCGGAGATG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

ORF Size:	1005 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001293697.2
RefSeq Size:	3275 bp
RefSeq ORF:	1008 bp
Locus ID:	10801
UniProt ID:	Q9UHD8
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome
MW:	38.5 kDa
Gene Summary:	This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2009]