

Product datasheet for **RG237556**

CLN3 (NM_001286105) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLN3 (NM_001286105) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CLN3
Synonyms:	BTN1; BTS; JNCL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237556 representing NM_001286105. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGACCCAGGCCAACGCGGATCCCCACAACAGCTCATCACGATTTGACTGCAACTCTGTCTCTA
CGGCTCCCCGGGTTCTCGTCAGTGGGATTTGTGCTGCTGGAAGCTTCGTCCTGGTTGCCTTTTCTCAT
TCTGTGGGACCAGCCTGTGTGGTGTGGTCTTCGCTAGCATCTCATCAGGCCTTGGGGAGGTCACCTTC
CTCTCCCTCACTGCCTTCTACCCAGGGCCGTGATCTCCTGGTGGTCCTCAGGGACTGGGGAGCTGGG
CTGCTGGGGCCCTGTCTACCTGGCCCTACCCAGGCCGCTCTCCCCTCAGCAGACCCTGCTGTCC
ATGCTGGGTATCCCTGCCCTGTGCTGGCCAGCTATTTCTTGTGCTCACATCTCCTGAGGCCCAGGAC
CCTGGAGGGGAAGAAGAAGCAGAGAGCGCAGCCCGCAGCCCCTATAAGAACCAGGCCCGGAGTCG
AAGCCAGGCTCCAGCTCCAGCCTCTCCCTTCGGGAAAGGTGGACAGTGTTCAAGGGTCTGCTGTGGTAC
ATTGTTCCCTTGGTCTAGTCTTACTTTGCCGAGTATTTTATTAAACCAGGGACTTTTTGAACCTCTTTT
TTCTGGAACACTTCCCTGAGTCACGCTCAGCAATACCGTGGTACCAGATGCTGTACCAGGCTGGCGTC
TTTGCCTCCCGCTTCTCTCCGCTGCTGTGCGATCCGTTTACCTGGCCCTGGCCCTGCTGCAGTGC
CTCAACCTGGTGTCTGCTGGCAGACGTGTGGTTCGGCTTTCTGCCAAGCATCTACCTCGTCTTCTCTG
ATCATTCTGTATGAGGGCTCCTGGGAGGCGCAGCCTACGTGAACACCTCCACAACATCGCCCTGGAG
ACCACTGATGAGCACCGGAGTTTGCAATGGCGGCCACCTGCATCTGACACACTGGGGATCTCCCTG
TCGGGGCTCCTGGCTTTGCCTCTGCATGACTTCTCTGCCAGCTCTCC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237556
 Blue=ORF Red=Cloning site Green=Tag(s)


MWTQAQRSPPTAHHDLTATLSLRLPRVLVSGICAAGSFVLVAFSHSVGTSLCGVVFASISSGLGEVTF
 LSLTAFYPRAVISWSSGTGGAGLLGALSYLGLTQAGLSPQQTLLSMLGIPALLLASYFLLLTSPEAQD
 PGEEEEAEASAARQPLIRTEAPESKPGSSSSLSLRERWTVFKGLLWYIVPLVVVYFAEYFVFNQGLFELLF
 FWNTSLSHAQQRWYQMLYQAGVFASRSSLRCCRIRFTWALALLQCLNLVFLLDVWVFGFLPSIYLVFL
 IILYEGLLGAAAYVNTFHNAIETSDEHREFAMAATCISDTLGLISLGLLALPLHDFLCQLS
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP
 SVIFTDKIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001286105

ORF Size:	1014 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001286105.2
RefSeq Size:	1717 bp
RefSeq ORF:	1017 bp
Locus ID:	1201
UniProt ID:	Q13286
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome
MW:	37.5 kDa
Gene Summary:	This gene encodes a protein that is involved in lysosomal function. Mutations in this, as well as other neuronal ceroid-lipofuscinosis (CLN) genes, cause neurodegenerative diseases commonly known as Batten disease or collectively known as neuronal ceroid lipofuscinoses (NCLs). Many alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]