

Product datasheet for **RG237567**

G protein alpha Inhibitor 2 (GNAI2) (NM_001282620) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	G protein alpha Inhibitor 2 (GNAI2) (NM_001282620) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	G protein alpha Inhibitor 2
Synonyms:	GIP; GNAI2B; H_LUCA15.1; H_LUCA16.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237567 representing NM_001282620. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGACGGAGGGTGTGAAGACGCTAGGCTGGACGAAGCAGAAAGCGGGTGTCACTGGGGACGTTCTGAG
GGTGTGGGGAGTCAGGGAAGAGCACCATCGTCAAGCAGATGAAGATCATCCACGAGGATGGCTACTCC
GAGGAGGAATGCCGCGAGTACCGGGCGGTTGTCTACAGCAACACCATCCAGTCCATCATGGCCATTGTC
AAAGCCATGGGCAACCTGCAGATCGACTTTGCCGACCCCTCCAGAGCGGACGACGCCAGGCAGCTATTT
GCACTGTCCTGCACCGCCGAGGAGCAAGGCGTGTCTCCCTGATGACCTGTCCGGCGTCATCCGGAGGCTC
TGGGCTGACCATGGTGTGCAGGCCTGCTTTGGCCGCTCAAGGGAATACCAGCTCAACGACTCAGCTGCC
TACTACCTGAACGACCTGGAGCGTATTGCACAGAGTGACTACATCCCCACACAGCAAGATGTGCTACGG
ACCCGCGTAAAGACCACGGGGATCGTGGAGACACACTTCACTTCAAGGACCTCACTTCAAGATGTTT
GATGTGGGTGGTCAGCGGTCTGAGCGGAAGAAGTGGATCCACTGCTTTGAGGGCGTCACAGCCATCATC
TTCTGCGTAGCCTTGAGCGCCTATGACTTGGTGTAGCTGAGGACGAGGAGATGAACCGCATGCATGAG
AGCATGAAGCTATTTCGATAGCATCTGCAACAACAAGTGGTTCACAGACACGTCCATCATCTCTTCCCTC
ACAAGAAGGACCTGTTTGAGGAGAAGATCACACACAGTCCCCTGACCATCTGCTCCCTGAGTACACA
GGGGCCAAACAAATATGATGAGGCAGCCAGCTACATCCAGAGTAAGTTTGAGGACCTGAATAAGCGCAAA
GACACCAAGGAGATCTACAGCACTTACAGTGCAGCCACCGACCAAGAACGTGCAGTTCGTGTTTGAC
GCCGTCACCGATGTCATCATCAAGAACAACCTGAAGGACTGCGGCCTCTTC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237567
 Blue=ORF Red=Cloning site Green=Tag(s)

MTEGVKTLGWTKQKGGCHWGRSEGAGESGKSTIVKQMKIIHEDGYSEEECRQYRAVVYSNTIQSIMAIV
 KAMGNLQIDFADPSRADDARQLFALSCTAEEQGVLPDDL SGVIRRLWADHGVQACFGRSREYQLNDSAA
 YYLNDLERIAQSDYIPTQQDVLRTVRKTTGIVETHFTFKDLHFKMFDVGGQRSERKKWIHCFEGVTAII
 FCVALSAYDLVLAEDEEMNRMHESMKLFDSICNNKWFDTDSIILFLNKKDLFEEKITHSPLTICFPEYT
 GANKYDEAASYIQSKFEDLNKRKDTKEIYTHFTCATDTKNVQVFVDAVTDVVIKNNLKDCGLF
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMKSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001282620

ORF Size: 1017 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


RefSeq: [NM_001282620.2](#)

RefSeq Size: 2467 bp

RefSeq ORF: 1020 bp

Locus ID:	2771
UniProt ID:	P04899
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance, Chemokine signaling pathway, Gap junction, Leukocyte transendothelial migration, Long-term depression, Melanogenesis, Progesterone-mediated oocyte maturation, Tight junction
MW:	39.2 kDa
Gene Summary:	The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]

Product images:



Circular map for RG237567