

Product datasheet for **RG237713**

CLN3 (NM_001286109) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLN3 (NM_001286109) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CLN3
Synonyms:	BTN1; BTS; JNCL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237713 representing NM_001286109. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGAGTGCCGCCACGACATCCTTAGCCACAAGAGGACATCGGGAAACCAGAGCCATGCTGTGCTC
CTGGCGGACATCCTCCCCACACTCGTCATCAAATTGTTGGCTCCTTTGGCCTTCACCTGCTGCCCTAC
AGCCCCGGGTTCTCGTCAGTGGGATTTGTGCTGCTGGAAGCTTCGTCTGGTTGCCTTTTCTATTCT
GTGGGGACCAGCCTGTGTGGTGTGGTCTTCGCTAGCATCTCATCAGGCCTTGGGGAGGTACCTTCCTC
TCCCTCACTGCCTTCTACCCAGGGCCGTGATCTCCTGGTGGTCTCAGGGACTGGGGAGCTGGGCTG
CTGGGGCCCTGTCTACCTGGGCTCACCCAGGCCGCTCTCCCTCAGCAGACCCTGTGTCCATG
CTGGGTATCCCTGCCCTGCTGCTGGCCAGCTATTTCTTGTGCTCACATCTCCTGAGGCCAGGACCCT
GGAGGGGAAGAAGAAGCAGAGAGCGCAGCCCGCAGCCCTCATAGAACCAGGCCCGGAGTCAAG
CCAGGCTCCAGCTCCAGCCTCTCCCTTCGGGAAAGGTGGACAGTGTCAAGGGTCTGCTGTGGTACATT
GTTCCCTTGGTCGTAGTTACTTTGCCGAGTATTTCAATTAACCAGGGACTTTTTGAACCTCTTTTTTC
TGGAACTTCCCTGAGTCACGCTCAGCAATACCGCTGGTACCAGATGCTGTACCAGGCTGGCGTCTTT
GCCTCCCGCTCTCTCCGCTGCTGTGCGATCCGTTTACCTGGGCCCTGGCCCTGCTGCAGTGCCTC
AACCTGGTGTTCCTGCTGGCAGACGTGTGGTTCGGCTTTCTGCCAAGCATCTACCTCGTCTTCTGATC
ATTCTGTATGAGGGCTCCTGGGAGGCGCAGCCTACGTGAACACCTCCACAACATCGCCCTGGAGACC
AGTGATGAGCACCGGGAGTTTGCAATGGCGGCCACCTGCATCTCTGACACACTGGGGATCTCCCTGTCC
GGGCTCTGGCTTTGCCTCTGCATGACTTCTCTGCCAGCTCTCC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237713
 Blue=ORF Red=Cloning site Green=Tag(s)


MLSAHDILSHKRTSGNQSHAVLLADILPTLVIKLLAPLGLHLLPYSRVLVSGICAAGSFVLVAFSHS
 VGTSLCGVVFASISSGLGEVTFLLSLTAFYPRAVISWSSGTGGAGLLGALSYLGLTQAGLSPQQLLSM
 LGIPALLLASYFLLLTSPEAQDPGGEAAESAARQPLIRTEAPESKPGSSSLSLRERWTVFKGLLWYI
 VPLVVVYFAEYFINQGLFELLFWNTLSHAQQYRWYQMLYQAGVFASSSLRCCRIRFTWALALLQCL
 NLVFLADVWFGFLPSIYLVFLIILYEGLLGGAAYVNTFHNIALETSEHREFAMAATCISDTLGISLS
 GLLALPLHDFLCQLS
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPGYPENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001286109
ORF Size:	1080 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001286109.2
RefSeq Size:	1791 bp
RefSeq ORF:	1083 bp
Locus ID:	1201
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome
MW:	39.6 kDa
Gene Summary:	This gene encodes a protein that is involved in lysosomal function. Mutations in this, as well as other neuronal ceroid-lipofuscinosis (CLN) genes, cause neurodegenerative diseases commonly known as Batten disease or collectively known as neuronal ceroid lipofuscinoses (NCLs). Many alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]