

## Product datasheet for **RG237752**

### CLCC1 (NM\_001278203) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CLCC1 (NM_001278203) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CLCC1
Synonyms:	MCLC; RP32
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237752 representing NM_001278203. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGTGTTCTTTGCTCCTTTGTGAATGTCTGTTGCTGGTAGCTGGTTATGCTCATGATGACTGG
ATTGACCCACAGACATGCTTAACATGATGCTGCTTCAGGAACAATGAGAAAATCTCAGGCAAAATAT
GGTATTTCCAGGGAAAAGGATGTCAGTCTGACTTGTCTGCTGATGAAATATCAGAATGTTATCAC
AACTTGATTCTTTAACTTATAAGATTGATGAGTGTGAAAAGAAAAGAGGGAAGACTATGAAAGTCAA
AGCAATCCTGTTTTAGGAGATACTTAAATAAGATTTAATTGAAGCTGGAAAGCTTGGACTTGCACCT
GCAGTTACATTCACCACATTTGTAAACGGAGCCATTGAAGCATATTGAAAAGGAAGCTGGGAATTTATT
AAAGCACTCATGAAGGAAATCCAGCGCTGCTTCATCTTCCAGTCTGATAATTATGGCATTAGCCATC
CTGAGTTTCTGCTATGGTCTGGAATAACAGTTTCAATGTGCTGAGACATATAGGCGGTCTGAGAGCGAA
CCTCCCCAGGCACTTCGGCCACGGGATAGAAGACGGCAGGAGGAAATTGATTATAGACCTGATGGTGGG
GCAGGTGATGCCGATTTCCATTATAGGGGCCAAATGGGCCCACTGAGCAAGGCCCTTATGCCAAAACG
TATGAGGGTAGAAGAGAGATTTGAGAGAGAGAGATGTTGACTTGAGATTTGAGACTGGCAACAAGAGC
CCTGAAGTGTCCGGCATTGATGTACCAGACGCAGAGGCACGAGAGCATCCCACGGTGGTACCCAGT
CATAAATCACCTGTTTTGGATACAAAGCCCAAGGAGACAGGTGGAATCCTGGGGGAAGGCACACCCGAAA
GAAAGCAGTACTGAAAGCAGCCAGTCGGCCAAGCCTGTCTGGCCAAGACACATCAGGGAATACAGAA
GGTTCACCCGCAGCGGAAAAGGCCAGCTCAAGTCTGAAGCCGCAGGCAGCCAGACCAAGGCAGCACA
TACAGCCCGCAAGAGGTGTGGCTGGACCAGTGGACAGGATCCGGTCTAGCAGCCCTGTGGC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

**Protein Sequence:** >Peptide sequence encoded by RG237752  
 Blue=ORF Red=Cloning site Green=Tag(s)


MLCSLLLCECLLLVAGYAHDDDWIDPTDMLNYDAASGTMRSQAKYGISGEKDVSPDLSCADEISECYH  
 KLDSLTYKIDECEKKKREDYESQSNPVFRRYLNKILIEAGKGLALAVTFTTFVTEPLKHIGKGTGEFI  
 KALMKEIPALLHLPVLIIMALAILSFYCYGAGKSVHVLRHIGGPESEPPQALRPRDRRRQEEIDYRPDGG  
 AGDADFHYRGQMGPTEQGPYAKTYEGRREILRERDVLRFQTGNKSPEVLRADFVDPDAEAREHPTVVPS  
 HKSPVLDTKPKETGGILGEGTPKESSTESSQSAKPVSGQDTSGNTEGSPA AEKAQLKSEAAGSPDQGST  
 YSPARGVAGPRGQDPVSSPCG  
 TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV  
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED  
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA  
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001278203
<b>ORF Size:</b>	1098 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>RefSeq:</b>	<a href="#">NM_001278203.1</a> , <a href="#">NP_001265132.1</a>
<b>RefSeq Size:</b>	4151 bp
<b>RefSeq ORF:</b>	1101 bp
<b>Locus ID:</b>	23155
<b>UniProt ID:</b>	<a href="#">Q96S66</a>
<b>Cytogenetics:</b>	1p13.3
<b>Protein Families:</b>	Ion Channels: Other, Transmembrane
<b>MW:</b>	40.3 kDa
<b>Gene Summary:</b>	Seems to act as a chloride ion channel.[UniProtKB/Swiss-Prot Function]