

Product datasheet for **RG237877**

LASS3 (CERS3) (NM_001290343) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LASS3 (CERS3) (NM_001290343) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LASS3
Synonyms:	ARCI9; LASS3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237877 representing NM_001290343. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTTGGACGTTTAAAGAATGGTTCTGGTTGGAAAGATTCTGGCTTCCTCCAACAATAAAGTGGTCA
GATCTTGAGGATCAGCATGGACTCGTCTTTGTAACCTTCTCATTTATACGTGACAATCCATATGCT
TTTCTCTGCTGATTATCAGGCGTGTATTTGAAAAATTTGTTGCTTCACCTCTAGCAAAATCATTGGC
ATTAAGAGACAGTTCGAAAGGTTACACCAAACTACTGTCTTAGAGAATTTTTTCAAACATCCACAAGG
CAACCATTGCAAACCTGATATTTATGGACTGGCAAAGAAGTGAACCTGACGGAGCGCCAGGTGAAAGA
TGTTTTAGGAGTCGGCGGAATCAAGAGAGGCCCTCCAGGCTGAAGAAATCCAGGAAGCTTCTGGAGA
TTTGCATTTTACTTAATGATCACTGTTGCTGGAATTGCGTTTCTTTATGATAAACCTTGGCTATATGAC
TTATGGGAGGTTTGAATGGCTATCCCAAACAGCCCCTGCTGCCATCCAGTACTGGTACTACATTTTA
GAAATGAGTTTTTATTGGTCTCTGTTATTTAGACTTGGCTTTGATGTCAAGAGAAAGGATTTTCTAGCT
CATATCATCCACCACCTGGCTGCTATTAGTCTGATGAGCTTCTCTTGGTGTGCTAATTATATTCGCAGT
GGGACCCTCGTGATGATTGTACACGATGTGGCTGACATTTGGCTGGAGTCTGCTAAGATGTTTTCTTAT
GCTGGATGGACGACACCTGTAACACCCTGTTTTTCATCTTCTCCACCATATTTTTCATCAGCCGCCCTC
ATTGTTTTTCTTCTGGATTTTATATTGCACGCTGATCTTGCCTATGTATCACCTCGAGCCTTTCTTT
TCATACATCTTCTCAACCTACAGCTCATGATCTTGCAGGTCCTCACCTTTACCTGGGTTATTACATC
TTGAAGATGCTCAACAGATGTATATTCATGAAGAGCATCCAGGATGTGAGGAGTGTGACGAGGATTAT
GAAGAGGAAGAGGAAGAGGAAGAAGAAGAGGCTACCAAAGGCAAAGAGATGGATTGTTTAAAGAACGGC
CTCGGGGCTGAGAGGCACCTCATTCCCAATGGCCAGCATGGCCAT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237877
 Blue=ORF Red=Cloning site Green=Tag(s)

MFWTFKEFWLERFWLPPTIKWSDLEDHDLVFKPSHLVYVIPYAFLLLIIRRVFEKFFVASPLAKSFG
 IKETVRKVTPNTVLENFFKHSTRQPLQTDIYGLAKKCNLTERQVERWFRSRNQERPSRLKKFQEACWR
 FAFYLMITVAGIAFLYDKPWLIDLWEVWNGYPKQPLLPSSQYWYIEMSFYWSLLFRLGFDVKKRDFLA
 HIIHHLAAISLMSFSWCANYIRSGTLVMIVHDVADIWLESAMFYSAGWTQTCNTLFFIFSTIFFISRL
 IVFPFWILYCTLILPMYHLEPFPSYIFLNLQLMILQVLHLYWGYIYLKMLNRCIFMKSIQDVRSDDEDY
 EEEEEEEEEATKGEMDCKNGLGAERHLIPNGQHG
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001290343

ORF Size: 1149 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


RefSeq: [NM_001290343.1](#), [NP_001277272.1](#)

RefSeq Size: 3913 bp

RefSeq ORF: 1152 bp
Locus ID: 204219
UniProt ID: [Q8IU89](#)
Cytogenetics: 15q26.3
Protein Families: Transcription Factors, Transmembrane
MW: 46.2 kDa

Gene Summary: This gene is a member of the ceramide synthase family of genes. The ceramide synthase enzymes regulate sphingolipid synthesis by catalyzing the formation of ceramides from sphingoid base and acyl-coA substrates. This family member is involved in the synthesis of ceramides with ultra-long-chain acyl moieties (ULC-Cers), important to the epidermis in its role in creating a protective barrier from the environment. The protein encoded by this gene has also been implicated in modification of the lipid structures required for spermatogenesis. Mutations in this gene have been associated with male fertility defects, and epidermal defects, including ichthyosis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2015]

Product images:



Circular map for RG237877