

Product datasheet for **RG237889**

Centaurin alpha 1 (ADAP1) (NM_001284308) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Centaurin alpha 1 (ADAP1) (NM_001284308) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Centaurin alpha 1
Synonyms:	CENTA1; GCS1L; p42IP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237889 representing NM_001284308. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGTTCCAGTTTGTTCAGCCGGGTGACTGCATTAACCTGCAAGAAGAAAATGGAAAGAATTTGAA
AAGATGCTGGGGTGTGCAGAGGAAGGGCACGCGAGTCTGGGAAGAGATCCCGACTGGGCCTCTACACT
CTGGGCGTCTTCATCTGCCTGAGCTGCTCGGGAATCCACCAGGAATATCCCCAGGTCAGCAAGGTGAAG
TCCGTCGGCTGGACGCCTGGGAGGAGGCCAAGTGGAGTTCATGGCCTCCACGGGAACGACGCCCGG
AGAGCCAGGTTT GAGTCCAAAGTACCCTCTTCTACTACCGGCCACGCCCTCCGACTGCCAGCTCCTT
CGAGAGCAGTGGATCCGGCCAAGTACGAGCGACAGGAGTTCATCTACCGGAGAAGCAGGAGCCCTAC
TCGGCAGGGTACCGT GAGGGTTTTCTCTGGAAGCGTGGCCGGGACAACGGGCAGTTTTT GAGCCGGAAG
TTTGTGCTGACAGAACGAGAGGGTGCTCTGAAGTATTTCAACAGAAATGATGCCAAGGAGCCCAAGGCC
GTGATGAAGATCGAGCACCTGAACGCCACCTTCCAGCCGGCCAAGATCGGCCACCCCCACGGCCTGCAG
GTCACCTACCTGAAGGACAACAGCACCCGTAACATCTTCATCTACCATGAGGACGGGAAGGAGATTGTG
GACTGGTTCAATGCACTCCGAGCTGCTCGCTTCCACTACCTGCAGGTGGCATTCCAGGGGCCAGCGAC
GCAGATCTGGTGCCAAAGCTCTCCAGGAACCTGAAGGAAGGCTACATGGAGAAGACGGGGCCCAAG
CAAACGGAAAGGCTTCCGGAAGCGCTGGTTCACCATGGATGACCGCAGGCTCATGTACTTCAAAGACCCC
CTGGACGCCTTCGCCGAGGGGAAGTCTTCATTGGCAGCAAGGAGTGGCTACACGGTGCATGGG
TTCCCGCGTCCACCCAGGGCCACCACTGGCCACATGGCATCACCATCGTCACGCCGACCGCAAGTTT
CTGTTTGCCTGCGAGACGGAGTCCGACCAGAGGAGTGGGTGGCGCCTTCCAGAAGCGGTGGACAGG
CCCATGCTGCCCCAGGAGTACGCAGTGGAGGCGCACTTCAAGCATAAACCT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237889
 Blue=ORF Red=Cloning site Green=Tag(s)

MFQVFVSRVYICNPARRKWEFEKMLGCAEEGHASLGRDPDWASYTLGVFICLSCSGIHRNIPQVSKVK
 SVRLDAWEEAQVEFMASHGNDAARARFESKVPFSFYRPTSDCQLLREQWIRAKYERQEFIYPEKQEPY
 SAGYREGFLWKRGRDNGQFLSRKFVLTREGALKYFNNDAKEPKAVMKIEHLNATFQPAKIGHPHGLQ
 VTYLKDNRNIFIYHEDGKEIVDFWNLRAARFHYLQVAFPGASDADLVKLSRNYLKEGYMEKTGPK
 QTEGFRKRWFTMDDRRLMYFKDPLDAFARGEVFIGSKESGYTVLHGFPPSTQGHHPHGITIVTPDRKF
 LFACETESDQREWVAAFQKAVDRPMLPQEYAVEAHFKHKP
 TRTRPLEMESDESLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPYSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNTAVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001284308

ORF Size: 1155 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NM_001284308.1](#), [NP_001271237.1](#)

RefSeq Size: 2195 bp

RefSeq ORF:	1158 bp
Locus ID:	11033
UniProt ID:	<u>O75689</u>
Cytogenetics:	7p22.3
Protein Families:	Druggable Genome
MW:	45.4 kDa
Gene Summary:	GTPase-activating protein for the ADP ribosylation factor family (Probable). Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG237889