

Product datasheet for **RG237913**

Cytochrome P450 3A5 (CYP3A5) (NM_001291829) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cytochrome P450 3A5 (CYP3A5) (NM_001291829) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CYP3A5
Synonyms:	CP35; CYP11A5; P450PCN3; PCN3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237913 representing NM_001291829. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAAAGTGCCATCTCTTTAGCTGAGGATGAAGAATGGAAGAGAATACGGTCATTGCTGTCTCCAACC
TTCACCAGCGGAAAACCAAGGAGATGTTCCCATCATTGCCAGTATGGAGATGTATTGGTGAGAAAAC
TTGAGGCGGGAAGCAGAGAAAAGCAAGCCTGTCACCTGAAAGACATCTTTGGGGCTACAGCATGGAT
GTGATTACTGGCACATCATTGGAGTGAACATCGACTCTCTCAACAATCCACAAGACCCCTTTGTGGAG
AGCACTAAGAAGTTCCTAAAATTTGGTTTCTTAGATCCATTATTTCTCTCAATAACTCTTTCCATTTC
CTTACCCAGTTTTTGAAGCATTAAATGTCTCTGTTTCCAAAAGATACCATAAATTTTTAAAGTAAA
TCTGTAACAGAAATGAAGAAAAGTCGCCTCAACGACAAAACAAAAGCACCGACTAGATTTCCCTCAGCTG
ATGATTGACTCCCAGAATTCGAAAGAAAAGTCCACAAAAGCTCTGTCTGATCTGGAGCTCGCAGCC
CAGTCAATAATCTTCATTTTTGCTGGCTATGAAACCACCAGCAGTGTCTTTCTTCACTTTATATGAA
CTGGCCACTCACCTGATGTCCAGCAGAAAAGTCAAAAAGGAGATTGATGCAGTTTTGCCAATAAGGCA
CCACCTACCTATGATGCCGTGGTACAGATGGAGTACCTTGACATGGTGGTGAATGAAACTCAGATTA
TTCCAGTTGCTATTAGACTTGAGAGGACTTGAAGAAAGATGTTGAAATCAATGGGGTATTCATTTCCC
AAAGGGTCAATGGTGGTATTCCAATTATGCTCTTACCATGACCCAAAGTACTGGACAGAGCCTTGAG
GAGTTCGGCCCTGAAAGTTTCAGTAAGAAGAAGGACAGCATAGATCCTTACATATACACACCCCTTTGGA
ACTGGACCCAGAAAACGATTGGCATGAGGTTTGTCTCATGAACATGAAACTTGCTCTAATCAGAGTC
CTTCAGAACTTCTCCTTCAAACCTGTAAAGAAAACACAGATCCCCTTGAATTAGACACGCAAGGACTT
CTTCAACCAGAAAACCCATTGTTCTAAAGGTGGATTCAAGAGATGGAACCCTAAGTGGAGAA
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237913
 Blue=ORF Red=Cloning site Green=Tag(s)


MKSAISLAEEDEEWKIRISLLSPTFTSGKLEMPFIIAQYGDVLRNLRREAEKGPVTLKIDIFGAYSMD
 VITGTSFGVNIIDSLNNPQDPFVESTKKFLKFGFLDPLFLSIIILFPFLTPVFEALNVSLFPKDTINFLSK
 SVNRMKKSRLNDKQKHRLDLFLQMLDSQNSKETESHKALSDLELAAQSIIFIFAGYETTSSVLSFTLYE
 LATHPDVQKQLQKEIDAVLPNKAPPTYDAVVQMEYLDMVVNETLRLFPVAIRLERTCKKDVVEINGVFIP
 KGSVMVVIPTYALHHDPKYWTEPEEFRPERFSKKKDSIDPYIYTPFGTGRNCIGMRFALMNMKLALIRV
 LQNF SFKPKETQIPLKLDTQGLLQPEKPIVLKVDSRDGTLSGE
TRTRPLEME SDE SGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYP SGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001291829
ORF Size:	1167 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001291829.2
RefSeq Size:	2255 bp
RefSeq ORF:	1170 bp
Locus ID:	1577
Cytogenetics:	7q22.1
Protein Families:	Druggable Genome, P450, Transmembrane
Protein Pathways:	Drug metabolism - cytochrome P450, Drug metabolism - other enzymes, Linoleic acid metabolism, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Retinol metabolism
MW:	44.8 kDa
Gene Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The encoded protein metabolizes drugs as well as the steroid hormones testosterone and progesterone. This gene is part of a cluster of cytochrome P450 genes on chromosome 7q21.1. Two pseudogenes of this gene have been identified within this cluster on chromosome 7. Expression of this gene is widely variable among populations, and a single nucleotide polymorphism that affects transcript splicing has been associated with susceptibility to hypertension. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]