

## Product datasheet for **RG237917**

### ABI2 (NM\_001282932) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ABI2 (NM_001282932) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ABI2
Synonyms:	ABI-2; ABI2B; AblBP3; AIP-1; AIP1; argBP1; argBP1A; argBPIB; SSH3BP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237917 representing NM_001282932. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTTTTAAAAGTTATTTACTTAAAATCTATACTTAGGCAAATTAGAGGCGTTGATCTTGAGTCG
ACTTTTGTGACCAAATTTGAAACAATTGCAGTTTGTGAGATTGAATGAGACAGTTGATATTCATAAAGAG
AAAGTTGCAAGAAGAGAAATTTGGTATTTTGTACTACCAATAAAAACACTTCAAGGACACATAAGATTATT
GCTCCAGCCAACCTTGAACGACCAGTTCGTTATATTAGAAAACCTATTGACTATAACAATTCTAGATGAT
ATTGGACATGGAGTAAAGGTGAGTACCCAGAACATGAAGATGGTGGGCTGCCGCTACAACACCTCCA
ACTCAGAAAGCCCTAGTCCCCCTATGTCAGGGAAAGGGACACTTGGGCGCACTCCCCCTATCGCACA
CTGGAGCCAGTGCCTCCAGTGGTACCAAATGATTACGTACCTAGCCCAACCGTAATATGGCTCCC
TCGCAGCAGAGCCCTGTGAGGACAGCTTCTGTGAATCAAAGAAATCGAACTTACAGCAGCAGTGGGAGT
AGTGGAGGGAGCCACCCAAGTAGTCGGAGCAGCAGTCGAGAGAACAGTGGAAAGTGGTAGTGTGGGGTT
CCTATTGCTGTTCTACTCCATCTCCTCCAGTGTCTTTCCAGGTCATCTGTACAGTTCTACAGCATG
AATAGGCTGCCTCTCGCCATACTCCCCAACAAATAGGGGGCTCGTTGCCCTATAGACGCCCTCCTTCC
ATTACTTCACAAACAAGCCTTCAAGATCAGATGAATGGAGGACCTTTTATAGCCAGAATCCAGTTTCA
GATACACCACCTCCACCGCCACCTGTGGAAGAACCAGTCTTTGATGAGTCTCCCCACCTCCTCCTCCT
CCAGAAGATTACGAAGAGGAGGAAGCTGTGGTTGAGTATAGTACCTTATGCTGAAGAGGACCCA
CCGTGGGCTCCACGTTCTTACTTGGAAAAGGTTGTGGCAATTTATGACTATACAAAAGACAAGGAAGAT
GAGCTGTCTTTTCAAGGAAGGAGCCATTATTTATGTCATCAAGAAGAATGACGATGGTTGGTATGAGGGA
GTTATGAATGGAGTACTGGGCTTTTTCTGGGAATTACGTTGAGTCTATCATGCATTATTCTGAG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

**Protein Sequence:** >Peptide sequence encoded by RG237917  
 Blue=ORF Red=Cloning site Green=Tag(s)


MFFKSYLLKISILRQIRGVDLESTFVTKFGNNCSLRLNETVDIHKEKVARREIGILTTNKNTSRTHKII  
 APANLERPVRYIRKPIDYITLDDIGHGKVVSTQNMKMGGLPRTTPPTQKPPSPMSSGKGLGRHSPYRT  
 LEPVRPPVVPNDYVPSPTRNMAPSQSPVRTASVNQRNRTYSSSSGSSGGSHPPSSRSSRENSGSGSVGV  
 PIAVPTSPSPSVFPGHPVQFYSMNRPASRHTPPTIGGSLPYRRPPSITSQTSLQNMNGGPFYSONPVS  
 DTTPPPPPVEEPVFDESPPPPPPEDEEEEEAAVVEYSDPYAEEDPPWAPRSYLEKVVAIYDYTKDKED  
 ELSFQEGAIYVIKKNDDGWYEGVMNGVTGLFPGNYVESIMHYSE  
**TRTRPLEMESDEGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV**  
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPEP  
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA  
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001282932
<b>ORF Size:</b>	1170 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>RefSeq:</b>	<a href="#">NM_001282932.2</a>
<b>RefSeq Size:</b>	6643 bp
<b>RefSeq ORF:</b>	1173 bp
<b>Locus ID:</b>	10152
<b>Cytogenetics:</b>	2q33.2
<b>Protein Pathways:</b>	Regulation of actin cytoskeleton
<b>MW:</b>	43.5 kDa
<b>Gene Summary:</b>	Regulator of actin cytoskeleton dynamics underlying cell motility and adhesion. Functions as a component of the WAVE complex, which activates actin nucleating machinery Arp2/3 to drive lamellipodia formation (PubMed:21107423). Acts as regulator and substrate of nonreceptor tyrosine kinases ABL1 and ABL2 involved in processes linked to cell growth and differentiation. Positively regulates ABL1-mediated phosphorylation of ENAH, which is required for proper polymerization of nucleated actin filaments at the leading edge (PubMed:7590236, PubMed:8649853, PubMed:10498863). Contributes to the regulation of actin assembly at the tips of neuron projections. In particular, controls dendritic spine morphogenesis and may promote dendritic spine specification toward large mushroom-type spines known as repositories of memory in the brain (By similarity). In hippocampal neurons, may mediate actin-dependent BDNF-NTRK2 early endocytic trafficking that triggers dendrite outgrowth (By similarity). Participates in ocular lens morphogenesis, likely by regulating lamellipodia-driven adherens junction formation at the epithelial cell-secondary lens fiber interface (By similarity). Also required for nascent adherens junction assembly in epithelial cells (PubMed:15572692).[UniProtKB/Swiss-Prot Function]