

Product datasheet for **RG237939**

TSPAN10 (NM_001290212) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TSPAN10 (NM_001290212) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TSPAN10
Synonyms:	OCSP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237939 representing NM_001290212. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCACAGACTAGCCAGGGCCGCGAGTCAGGACCAGCCAGCAGCCAGCCACAGAGGAGCCAGCGAAC
CTCTCCCGGCGCCTGTGCTGGGGCTTTCTGTTCCAGCGTCAAGGATGGAGGAGGGGAGAGGAGCCCC
TTACTGTCCCAGGAACTGCAGGCCAGAAGCCCTCTCTGTGCACAGGCCACCCACCTCAGGCTGCCTA
GGTCCAGTGCCAGGGAGGACCAGGCGGAGGCCTGGGGCTGCAGCTGCTGTCCCGGAGACCAAGCAC
CAGGCCTTGAGTGGCACCCEAAGAAAGGACCAGCCCTTCCCTCTCCCAGGGAGCAGCTGCGTCAAG
TATCTGATCTTCTCTCCAACCTCCCTTCTCCCTGCTGGGGCTGCTGGCCCTGGCCATCGGGCTCTGG
GGCCTGGCTGTCAAGGGTCTCTGGGAAGTGATCTGGGGGGGCCCTGCCCGCAGACCCCATGCTGGGG
CTGGCACTGGGAGGGCTGGTGGTCAGCGCAGTGAGCCTGGCTGGCTACCTGGGCGCCCTCTGTGAGAAC
ACCTGCCTGTTACGTGGCTTCTCTGGGGCATCCTTGCCCTTCTGGTGTGAGGCCGTGGCGGGGGCC
CTGGTGGTGGCCCTCTGGGGCCCGCTGCAAGACAGCCTGGAGCACACCCTGCGTGTGGCCATCGCCAC
TACCAGGACGACCAGACCTGCGCTTCTCCTCGACCAAGTCCAGCTCGGGCTGAGGTGCTGCGGAGCT
GCCTCTACCAGGACTGGCAGCAGAACCTGTACTTTAACTGCAGCTCCCCGGGTGCAGGCCTGCAGC
CTTCCCGCCTCTGCTGCATCGACCCCGCAAGATGGAGCCTCTGTCAACGACCAGTGGCGCTTCGGG
GTCCTGCGCCTGGATGCGGACGACGCTCAGAGAGTGGTGTACCTGGAGGGCTGCGGCCGCGCCTCCGG
CGGTGGCTGCGCGCAACCTGGCTGCCTCGGGCGGCTACGCAATCGCGGTGGTGTGCTGCAGGGCGCG
GAGCTCTGCTGGCCGCCGGCTACTCGGGCCCTCGCTGCCCGCAGGGGGCGGCGTACGGCCCGGA
GCGCGGGGGAGGACCGCTGGCCCCAGAGCCAGCCCGCGCCCGCCGCTGCCAAGCCCGCC
CGGGC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```




[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237939
 Blue=ORF Red=Cloning site Green=Tag(s)


MHRLAQRSDQPSSPATEEPANLSRRLCWGLSVPASRMEEGERSPLLSQETAGQKPLSVHRPPTSGCL
 GPVPREDQAEAWGCSCCPPETKHQALSGTPKKGPAAPSLSPGSSCVKYLIFLSNPFLLGLLALAIGLW
 GLAVKGSLSGDLGGPLPADPMLGLALGGLVVSAYSLAGYLGALCENTCLLRGFSGGILAFVLVLEAVAGA
 LVVALWGPLQDSLEHTLRVAIAHYQDDPDLRFLLDQVQLGRCCGAASYDQWQNLVFNCSPPGVQACS
 LPASCCIDPREDGASVNDQCGFVLRRLDADAAQRVVYLEGCGPPLRRWLRANLAASGGYAI AVVLLQGA
 ELLLAARLLGALAARRGAAYGPGARGEDRAGPQSPSPGAPPAAKPARG
TRTRPLEMEDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPSTYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001290212
ORF Size:	1179 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001290212.1 , NP_001277141.1
RefSeq Size:	1893 bp
RefSeq ORF:	1068 bp
Locus ID:	83882
UniProt ID:	Q9H1Z9
Cytogenetics:	17q25.3
Protein Families:	Transmembrane
MW:	41.3 kDa
Gene Summary:	Regulates maturation of the transmembrane metalloprotease ADAM10.[UniProtKB/Swiss-Prot Function]