

Product datasheet for **RG237980**

ZNF302 (NM_001289182) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF302 (NM_001289182) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF302
Synonyms:	HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG237980 representing NM_001289182. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCT
GCTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGCAGGTCTT
TCCGTAACCAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGAGCCCTGGATGATGGAGAAAAA
CTGTCAAAGATTGGGAATCAAGATGGGAAAACAAGGAATTATCAACAAAGAAGGATATTTATGATGAA
GATTCACCCCAACAGTAACAATGGAAAAGTTGTAACAACAAGTTATGAATTTCAAATCTAATAAG
AATTTGGAATATACAGAATGCGACACATTTAGAAGCACCTTTCATTCAAAGTCTACTCTTTCTGAACCA
CAAAACAATTCTGCTGAAGGGAATTCACACAAATATGATATATTAAGAAGAATTTATCAAAAAAGTCA
GTTATAAAAAGTGAGAGAATAAATGGTGGAAAAGAACTTTTAAATCTAATAAAAAGTGGGCGAGCCTTC
AACCAGAGCAAATCTCTTACCCTTCCCCAGACTTGAATAGAGAGAAAATCTATACATGCAGTGAATGT
GGGAAAGCCTTTGGCAAACAGTCAATCCTCAGTCGCCACTGGAGAATTCATACAGGAGAGAAGCCCTAT
GAATGTCGTGAATGTGGGAAGACTTTTAGCCATGGTTCATCCCTTACACGACATCAGATAAGCCATAGT
GGAGAGAAACCTTACAAATGCATTGAATGTGGGAAGGCCTTTAGCCATGGCTCATCACTTACTAACCAT
CAGAGCACTCACACGGGAGAGAAAACCGTATGAATGTATGAACTGTGGAAGTCTTTTAGTCGTGTGTCC
CTTCTCATTGAGCATCTAAGAATTCATACGCAAGAAAACGCTATGAGTGTGATATGTGGAAAGGCC
TTTATTAGTTTCGTCTCTCATTCAACATCAGAAAAGCCATACTGGAGAGAAGCCTTATGAATGTAGA
GAATGTGGAAAAGCTTTCTGCTGTAGCTCACACCTTACTCAACATCAAAGAATTCACAGTATGAAGAAA
AAATATGAATGCAACAAATGTCTCAAGGCTTTTAGTAGCTTCTCATTTCTGTTCAACATCAGAGTATT
CATACTGAAGAAAACCGTTTGAAGTT
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
```



[View online »](#)

Protein Sequence: >Peptide sequence encoded by RG237980
 Blue=ORF Red=Cloning site Green=Tag(s)


MSQVTFSDVAIDFSHEEWACLDSAQRDL YKDV MVQNYENLVSVAGLSVTKPYVIMLLEDGKEPMMMEKK
 LSKDWESRWENKELSTKKDIYDEDSPQPVMEKVVKQSYEFSNSNKNLEYTECDTFRSTFHSKSTLSEP
 QNNSAEGNSHKYDILKKNLSKKSVIKSERINGGKLLNSNKSAAFNQSKSLTLPQTCNREKIYTCSEC
 GKAFGKQSILSRHWRIHTGEKPYECRECGKTF SHGSSL TRHQISHSGEKPYKCIECGKAF SHGSSL TNH
 QSTHTGEKPYECMNCGKFSRVSLLIQHLRIHTQEKRYECRICGKAF IHSSSLIHHQKSHTEKPYECR
 ECGKAFCCSSHLTQHQRISHMKKKYECNKCLKVFSFSLVQHQSIIHTEKPFVEV
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
 MGYGFYHFGTYPGYNPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
 SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFLRDGGYSSVVD SHMHFKSAIHPSILQNGGPMFA
 FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001289182
ORF Size:	1200 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001289182.1 , NP_001276111.1
RefSeq Size:	2938 bp
RefSeq ORF:	1203 bp
Locus ID:	55900
UniProt ID:	Q9NR11
Cytogenetics:	19q13.11
Protein Families:	Transcription Factors
MW:	46.5 kDa
Gene Summary:	This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]