

Product datasheet for **RG238956**

PRMT5 (NM_001282953) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRMT5 (NM_001282953) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PRMT5
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG238956 representing NM_001282953.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGGGGTCCGAACTCGGGGACGGAGAAGGGCAGACTAGTCATCCCGGAGAAGCAGGGGTTTGATTC
CTCTGCATGCCTGTCTCCATCCCGCTTCAAGAGGGAGTTTCATTCAGGAACCTGCTAAGAATCGGCC
GGTCCCCAGACAGATCAGACCTACTGCTGTAGGAAAGGGCTTCTCTGCTGCCCTTAATCAGGAAGAT
AACACCAACCTGGCCAGAGTTTTGACCAACCACATCCACACTGGCCATCACTCTTCCATGTTCTGGATG
CGGGTACCCTTGGTGGCACCAGAGGACCTGAGAGATGATAAATTGAGAAATGCACCAACTACACACACA
GAGGAGTACAGTGGGAGGAGAAAACGTGGATGTGGTGGCACAACCTCCGGACTTGTGTACTATAGT
AAGAGGATTGCAGTGGCTCTTGAATTGGGGCTGACCTCCCATCTAATCATGTCATTGATCGCTGGCTT
GGGGAGCCCATCAAAGCAGCCATCTCCCACTAGCATTCTCTGACCAATAAGAAGGGATTCTCTGTT
CTTTCTAAGATGCACCAGAGGCTCATCTCCGGCTCCTCAAGTTGGAGGTGCAGTTCATCATCACAGGC
ACCAACCACCACTCAGAGAAGGAGTTCTGCTCCTACCTCCAATACCTGGAATACTTAAGCCAGAACCGT
CCTCCACCTAATGCCTATGAACTCTTGGCAAGGGCTATGAAGACTATCTGCAGTCCCCGCTTCAGCCA
CTGATGGACAATCTGGAATCTCAGACATATGAAGTGTGAAAAGGACCCATCAATACTCTCAGTAC
CAGCAGGCCATCTATAAATGTCTGCTAGACCGAGTACCAGAAGAGGAGAAGGATACCAATGTCCAGGTA
CTGATGGTGTGGGAGCAGGACGGGACCCCTGGTGAACGCTTCCTGCGGGCAGCCAAGCAGGCCGAC
CGGCGGATAAAGCTGTATGCTGTGGAGAAAAACCCAAATGCCGTGGTACGCTAGAGAATGGCAGTTT
GAAGAAATGGGGAAGCCAAGTGACCGTAGTCTCATCAGACATGAGGGAATGGGTGGCTCCAGAGAAAGCA
GACATCATTGTCAGTGAAGTCTTGGGCTCATTGCTGACAATGAATTGTCGCCTGAGTGCCTGGATGGA
GCCCAGCACTTCTAAAAGATGATGGTGTGAGCATCCCCGGGAGTACACTTCTTCTGGCTCCCATC
TCTTCTCAAGCTGTACAATGAGGTCGAGCCTGTAGGAGAAGGACCGTACCCCTGAGGCCAGTTT
GAGATGCCTTATGTGGTACGGCTGCACAACCTCCACCAGCTCTCTGCACCCAGCCCTGTTTCACCTTC
AGCCATCCCAACAGAGATCCTATGATTGACAACAACCGCTATTGCACCTTGAATTTCTGTGGAGGTG
AACACAGTACTACATGGCTTTGCCGGCTACTTTGAGACTGTGCTTATCAGGACATCACTCTGAGTATC
CGTCCAGAGACTCACTCTCTGGGATGTTCTCATGGTTTCCCATCCTCTTCCCTATTAAGCAGCCATA
ACGGTACGTGAAGCCAAACCATCTGTGTGCGTTTCTGGCGATGCAGCAATCCAAGAAGGTGTGGTAT
GAGTGGGCTGTGACAGCACCAGTCTGTTCTGCTATTCATAACCCACAGGCCGCTCATATACCATTGGC
CTC
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAAAC
```

Protein Sequence:


>Peptide sequence encoded by RG238956
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MRGPNSGTEKGRLLVIPEKQGFDFLCMPVHFHPRFKREFIQEPAKNRPGPQTRSDLLL SGRFLLPLNQED
NTNLRVLTNHIHTGHSSMFWMRVPLVAPEDLRDDIIENAPTHTEEYSGEEKTWMWWHNFRTLCDYS
KRIAVALEIGADLPSNHVIDRWLGEPIKAAILPTSI FL TNKKGFPVL SKMHQRLIFRLLKLEVQFIITG
TNHHSEKEFCSYLQYLEYLSQNRPPNAYELFAKGYEDYLSPLQLMDNLESQTYEVFEKDPKIKYSQY
QQAIYKCLLDRVPEEEKDTNVQVLMVLGAGRGPLVNASLRAAKQADRRIKLYAVEKNPNNAVVTLENWQF
EEWGSQVTVVSSDMREWVAPEKADIIIVSELLGSFADNELSPECLDGAQHFLKDDGVSIPGEYTSFLAPI
SSSKLYNEVRACREKDRDPEAQFEMPYVRLHNFHQLSAPQPCFTFSHPNRDPMIDNNRYCTLEFPVEV
NTVLHGFAGYFETVLYQDITLSIRPETHSPGMFSWFPILFPIKQPI TVREGQICVRFWRCSNSKKVWY
EWAVTAPVCSAIHNPTGRSYTIGL
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIIRSNAIVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSHMFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282953

ORF Size: 1728 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001282953.2
RefSeq Size:	2494 bp
RefSeq ORF:	1731 bp
Locus ID:	10419
UniProt ID:	O14744
Cytogenetics:	14q11.2
Protein Families:	Stem cell - Pluripotency
MW:	66.3 kDa
Gene Summary:	<p>This gene encodes an enzyme that belongs to the methyltransferase family. The encoded protein catalyzes the transfer of methyl groups to the amino acid arginine, in target proteins that include histones, transcriptional elongation factors and the tumor suppressor p53. This gene plays a role in several cellular processes, including transcriptional regulation, and the assembly of small nuclear ribonucleoproteins. A pseudogene of this gene has been defined on chromosome 4. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]</p>