

## Product datasheet for **RG239140**

### **BMAL1 (ARNTL) (NM\_001297722) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BMAL1 (ARNTL) (NM_001297722) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ARNTL
Synonyms:	bHLHe5; BMAL1; BMAL1c; JAP3; MOP3; PASD3; TIC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG239140 representing NM\_001297722.  
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCAGACCAGAGAATGGACATTTCTCAACCATCAGTGATTTTCATGTCCCGGGCCCCACCGACCTG
CTTCCAGCTCTCTTGGTACCACTGGTGTGGATTGCAACCGCAAACGGAAAGGCAGCTCCACTGACTAC
CAAGAAAGCATGGACACAGACAAAAGATGACCCTCATGGAAGGTTAGAATATACAGAACCCAAGGAAGG
ATAAAAAATGCAAGGGAAGCTCACAGTCAGATTGAAAAGCGGCGTCGGGATAAAAATGAACAGTTTATA
GATGAATTGGCTTCTTGGTACCAACATGCAACGCAATGTCCAGGAAATTAGATAAACTTACTGTGCTA
AGGATGGCTGTTCCAGCACATGAAAACATTAAGAGGTGCCACCAATCCATACACAGAAGCAAACACTACAAA
CCAACCTTTCTATCAGACGATGAATTGAAACACCTCATTCTCAGGGCAGCAGATGGATTTTTGTTTGTG
GTAGGATGTGACCGAGGGAAGATACTCTTGTCTCAGAGTCTGTCTTCAAGATCCTCAACTACAGCCAG
AATGATCTGATTGGTCAGAGTTTGTGTTGACTACCTGCATCCTAAAGATATTGCCAAAGTCAAGGAGCAG
CTCTCCTCTCTGACACCGCACCCCGGGAGCGGCTCATAGATGCAAAAAGTGGACTTCCAGTTAAAAACA
GATATAACCCCTGGGCCATCTCGATTATGTTCTGGAGCAGCAGCTTCTTCTTCTGTAGGATGAAGTGT
AACAGGCCTTCAGTAAAGGTTGAAGACAAGGACTTCCCCTCTACCTGCTCAAAGAAAAAAGCAGATCGA
AAAAGCTTCTGCACAATCCACAGCACAGGCTATTTGAAAAGTGGCCACCCACAAAGATGGGGCTGGAT
GAAGACAACGAACCAGACAATGAGGGGTGAACCTCAGCTGCCTCGTCGCAATTGGACGACTGCATTCT
CATGTAGTTCCACAACCAAGTGAACGGGAAATCAGGGTGAATCTATGGAATATGTTTCTCGGCACGGC
ATAGATGGAAAGTTTGTGTTTGTAGACCAGAGGGCAACAGCTATTTGGCATATTTACCACAAGAAGCTT
CTAGGCACATCGTGTATGAATATTTCCACCAAGATGACATAGGACATCTTGAGAATGTGCATAGGCCAA
GTTTTACAGACGAGAGAAAAAATTACAATAATTGCTATAAAATTTAAAAATCAAAGATGCTCTTTTATC
ACACTACGGAGTCGATGGTTTCAAGTTTCAATGAACCTTGGACCAAGGAAGTAGAATATATTGTCTCACT
AACACTGTTGTTTAGCCAACGTCCTGGAAGGCGGGACCCAACCTTCCCACAGCTCACAGCATCCCCC
CACAGCATGGACAGCATGCTGCCCTCTGGAGAAGGTGGCCCAAAGAGGACCCACCCCACTGTTCCAGGG
ATTCCAGGGGAACCCGGGCTGGGGCAGGAAAAATAGGCCGAATGATTGCTGAGGAAATCATGGAATC
CACAGGATAAGAGGGTCACTCGCCTTCTAGCTGTGGCTCCAGCCATTGAACATCACGAGTACGCCTCCC
CCTGATGCCTCTTCTCCAGGAGCAAGAAGATTTAAATGGAGGGACTCCAGACATTCTTCCAGTGGC
CTACTATCAGGCCAGGCTCAGGAGAACCAGGTTATCCATATTCTGATAGTTCTTCTATTCTTGGTGAG
AACCCCCACATAGGTATAGACATGATTGACAACGACCAAGGATCAAGTAGTCCCAGTAATGATGAGGCA
GCAATGGCTGTCATCATGAGCCTCTTGAAGCAGATGCTGGACTGGTGGCCCTGTTGACTTTAGTGAC
TTGCCATGGCCGCTG
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAAAC
```

**Protein Sequence:**

>Peptide sequence encoded by RG239140  
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MADQRMDISSTISDFMSPGPTDLLSSSLGTSGVDCNRKRKGSSTDYQESMDTDKDDPHGRLEYTEHQGR
IKNAREAHSQIEKRRRDKMNSFIDELASLVPTCNAMSRKLDKLTVLRMAVQHMKTLRGATNPYTEANYK
PTFLSDDELKHLILRAADGLFVVGCDRGKILFVSESVFKILNYSQNDLIGQSLFDYLHPKDIKVKKEQ
LSSSDTAPRERLIDAKTGLPVKTDITPGPSRLCSGARRSFFCRMKNRPSVKVEDKDFPSTCSKKKADR
KSFCTIHSTGYLKSPPKTKMGLDEDNEPDNEGCNL SCLVAIGRLHSHVVPQPVNGEIRVKSMEYVSRHA
IDGKFFVVDQRATAILAYLPQELLGTSCYEYFHQDDIGHLAECHRQVLQTREKITTNCYKFKIKDGSFI
TLRSRWF S F M N P W T K E V E Y I V S T N T V V L A N V L E G G D P T F P Q L T A S P H S M D S M L P S G E G G P K R T H P T V P G
I P G G T R A G A G K I G R M I A E E I M E I H R I R G S S P S S G S S P L N I T S T P P P D A S P G G K I L N G T P D I P S S G
L L S G Q A Q E N P G Y P Y S D S S I L G E N P H I G I D M I D N D Q G S S P S N D E A A M A V I M S L L E A D A G L G P V D F S D
L P W P L
TRTRPLEMESDESGLPAMEIECRITGTLNGVEFELVGGEGTPEQGRMTNKMSTKGALTFSPYLLSHV
MGYGFYHFGTYPSGYENPFLHAINNGGYNTRIEKYEDGGVLHVSFSYRYEAGRVIGDFKVMGTGFPED
SVIFTDKIIRSNATVEHLHPMGDNDLDGSFTRTFSLRDGGYSSVVDSSHMFKSAIHPSILQNGGPMFA
FRRVEEDHSNTELGIVEYQHAFKTPDADAGEERV
```


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_001297722

ORF Size: 1878 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>RefSeq:</b>	<a href="#">NM_001297722.1</a> , <a href="#">NP_001284651.1</a>
<b>RefSeq Size:</b>	2811 bp
<b>RefSeq ORF:</b>	1881 bp
<b>Locus ID:</b>	406
<b>UniProt ID:</b>	<a href="#">O00327</a>
<b>Cytogenetics:</b>	11p15.3
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>Protein Pathways:</b>	Circadian rhythm - mammal
<b>MW:</b>	69.2 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a basic helix-loop-helix protein that forms a heterodimer with CLOCK. This heterodimer binds E-box enhancer elements upstream of Period (PER1, PER2, PER3) and Cryptochrome (CRY1, CRY2) genes and activates transcription of these genes. PER and CRY proteins heterodimerize and repress their own transcription by interacting in a feedback loop with CLOCK/ARNTL complexes. Defects in this gene have been linked to infertility, problems with gluconeogenesis and lipogenesis, and altered sleep patterns. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]