

Product datasheet for **RG239865**

POLR1B (NM_001282779) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POLR1B (NM_001282779) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	POLR1B
Synonyms:	A135; RPA2; RPA135; Rpo1-2; TCS4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG239865 representing NM_001282779.
 Blue=ORF Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGTGAATCCAAGCTTTGCAACTTACGTAACCTTCCCCACAAGCCCTCATTGAGCACCATGAGGAG
CGAGAGGAAATGGGGGGCTATTTTATAATCAATGGCATTGAAAAAGTCATCCGAATGTTGATTATGCCT
CGGAGAAAATTTCCATTGCAATGATAAGACCAAATGGAAAACAGAGGCCTGGTTATACTCAGTAT
GGAGTTTCAATGCACTGTGTGAGGGAAGAACATCCGCTGTCAATATGAACCTCCACTACTTGAAAAAT
GGCACTGTTATGTTGAACTTTATTTACCAAAAGAAGTGTCTTTCTTCTTTGGGATTTGCACTTAAG
GCCTTGTGAGCTTTTCTGATTATCAGATCTTTCAGGAGCTCATCAAAGGAAAAGAGGATGATTCTTTC
CTTAGGAACTCTGTTTCTCAGATGTTAAGGATTGTAATGGAAGAGGGTGTTCGACACAAAAACAGGTC
CTTAACACTAGGTGAATGCTTCAGAGTAAACTCAATGTTCTGACTGGTACCCAAATGAGCAAGCT
GCGGAGTTCCTGTTAACCAGTGCATCTGTATCCACTTGAATCCAATACTGAAAAGTTTTATATGCTT
TGTCTCATGACCGAAAGCTCTTTGCTTACGCAAAGGAGAGTGCATGGAGGACAATCCTGATAGTTTG
GTGAACAGGAAGTCTCACACCGGTCAGCTCTTCTTATGTTCTGAAGGAAAAACTGGAAGTTGG
TTAGTGTCTATTAATAAGCTTTTGATAAGAAGGCTCAGAAGACCAAGTGTTCATGAACACTGACAAT
TTGATGAGGATTTTACAATGGGCATAGACCTTACAAAACCAATTTGAATACCTTTTTGCTACTGGGAAT
CTGCGTTCTAAAACAGGTCTTGGCCTCCTACAAGATTCTGGACTTGTGTTGTGGCTGACAAGCTGAAC
TTCATACGCTACCTCTCCATTTCCGCTGCGTGCACAGAGGGGCTGATTTTGCCAAGATGAGGACCACC
ACAGTACGCAGGTGCTGCCAGAGTCTGGGGCTTCTTTGTCCCGTGCATACCCAGACGGGGAGCCC
TGTGGCCTGATGAACCACCTAACTGCCGTATGTGAGGTTGTACACAGTGTGTATACGGCATCTATT
CCAGCTTTACTGTCAACTTGGGGTCACTCCCATTGATGGAGCTCCCCACCGATCATACAGTGAATGC
TACCCTGTCTGCTGGACGGTGTCTAGTTGGCTGGGTGGATAAAGGATCTTGTCCAGGCATCGCAGAT
TCTCTTCTGTCATTTAAGGTGTTGAGAGAGAAAAGAAATCCTCCCTGGATGGAAGTGGTCTTTATACCC
ATGACAGGAAAACCAAGTCTGTACCCAGGATTGTTCTTTTTACCACTCCTTGTAGACTGGTACGGCCT
GTGCAGAACTTAGCATTGGGCAAAGAAGAGCTAATTGGAAGTATGGAACAGATCTTCATGAATGTCGCT
ATCTTTGAGGATGAAGTTTTGCTGGAGTTACCACACACCAGGAACTCTTCCACACAGCCTGCTGAGT
GTGATTGCCAACTTCATCCCTTCTCTGATCACAACCAGAGTCCACGGAACATGTACCAATGCCAGATG
GGTAAGCAAATATGGGCTTTCCACTTCTCACTTATCAAGACCGATCGGATAACAACTGTATCGTCTT
CAGACTCCTCAGAGTCCCTTGGTGAACCTCCATGTATGATTATTATGACATGGATAACTATCCAATT
GGGACCAATGCCATCGTTGCTGTGATTTCTTACTGGCTATGATATGGAAGATGCCATGATTGTGAAT
AAGGCCTCTTGGGAACGAGGCTTTGCCATGGAAGTGTCTACAAGTCTGAGTTTATAGACCTCTCTGAA
AAAATTAACAAGGAGATAGTAGCTGGTGTGGCATCAAACCTGGTGACCCACGCTTCTGCAGAAG
TTAGATGACGATGGATTGCCGTTTATAGGAGCAAACTGCAGTACGGAGATCCGTATTACAGCTACCTC
AACCTCAACACCGGGGAAAGTTTTGTGATGTAATAAGAGTAAAGAAAATGTGTTGTGGATAACATC
AAAGTGTGAGTAAATGACTGGGAGTGGAAAATCAAGTGTGTTGCATCACTATGAGAGTGCCTCGG
AACCAACTATCGGAGATAAATTTGCCAGTCCCATGGGCAGAGGGCATTAAAGCAGATTGTGGCCG
GCTGAGACATGCCTTTTACTGAGAGTGGGATGGTCCAGACATCTGTTCATCCCATGGTTTTCCA
TCCCGCATGACCATTGGGATGTTAATTGAGAGTATGGCCGGGAAGTCTGCAGCTTGTGATGGTCTCTGC
CATGATGTACACCTTTCATCTTCTCAGAGGAGAACTCGGCCTTAGAATACTTTGGTGAAGTGTAAAG
GCTGCTGGCTACAATTTCTATGGCACCGAGAGTTATATAGTGGCATCAGTGGGTAGAAGTGAAGCA
GACATCTTCATAGGAGTGGTTTATATCAGCGCTTACGCCATATGGTCTCAGACAAATTTCAAGTAAAG
ACAAGTGGAGCCCGAGACAGAGTACCAACCAGCCTATTGGGGGAAGAAATGTCCAGGGTGGAAATCCGT
TTTGGGGAGATGGAACGGGATGCGCTTTAGCTCATGGTACATCTTTCTCCTTCATGACCGCTCTTC
AACTGCTCAGATCGTTCGGTAGCCCATGTGTGTGAAGTGTGGCAGTTTACTCTCCTACTGTTGGAG
AAGCCACCCCTTCTTGGTCTGCCATGCGCAACAGAAAATACAAGTACTCTGTGTAGTGCAGTGCAG
ACTATCGATACTGTTTCTGTGCCTTATGTTTTTCGGTATTTTGTAGCTGAACTGGCAGCTATGAACATC
AAAGTGAAGTGGATGTTGTT
AGCGGACCGACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAAAC
  
```



Protein Sequence: >Peptide sequence encoded by RG239865
 Blue=ORF Red=Cloning site Green=Tag(s)

MVKS^{Blue}KL^{Blue}CN^{Blue}LR^{Blue}NL^{Blue}PP^{Blue}QAL^{Blue}IEH^{Blue}HEE^{Blue}AE^{Blue}EM^{Blue}GG^{Blue}YF^{Blue}I^{Blue}ING^{Blue}IE^{Blue}KV^{Blue}IR^{Blue}ML^{Blue}IM^{Blue}PR^{Blue}RF^{Blue}PI^{Blue}AM^{Blue}IR^{Blue}PK^{Blue}WK^{Blue}TR^{Blue}GP^{Blue}Y^{Blue}TQ^{Blue}
 G^{Blue}VS^{Blue}M^{Blue}HC^{Blue}V^{Blue}REE^{Blue}HS^{Blue}AV^{Blue}N^{Blue}M^{Blue}NL^{Blue}HY^{Blue}LE^{Blue}NG^{Blue}TV^{Blue}ML^{Blue}NF^{Blue}I^{Blue}Y^{Blue}R^{Blue}KE^{Blue}LF^{Blue}FL^{Blue}PL^{Blue}GF^{Blue}AL^{Blue}KAL^{Blue}V^{Blue}SF^{Blue}SD^{Blue}YQ^{Blue}IF^{Blue}QEL^{Blue}IK^{Blue}GK^{Blue}ED^{Blue}DS^{Blue}F^{Blue}
 LR^{Blue}NS^{Blue}V^{Blue}S^{Blue}Q^{Blue}ML^{Blue}R^{Blue}IV^{Blue}ME^{Blue}EG^{Blue}CS^{Blue}TQ^{Blue}KQ^{Blue}VL^{Blue}NYL^{Blue}GE^{Blue}CF^{Blue}RV^{Blue}KL^{Blue}N^{Blue}VP^{Blue}DW^{Blue}YP^{Blue}NE^{Blue}QAA^{Blue}E^{Blue}FL^{Blue}N^{Blue}QC^{Blue}IC^{Blue}I^{Blue}HL^{Blue}KS^{Blue}NT^{Blue}E^{Blue}K^{Blue}F^{Blue}Y^{Blue}ML^{Blue}
 CL^{Blue}M^{Blue}TR^{Blue}KL^{Blue}F^{Blue}AL^{Blue}AK^{Blue}GE^{Blue}CM^{Blue}ED^{Blue}NP^{Blue}DSL^{Blue}V^{Blue}NQ^{Blue}EVL^{Blue}TP^{Blue}GQL^{Blue}FL^{Blue}M^{Blue}FL^{Blue}KE^{Blue}KE^{Blue}LE^{Blue}GL^{Blue}W^{Blue}VS^{Blue}IK^{Blue}IA^{Blue}FD^{Blue}KA^{Blue}QK^{Blue}TS^{Blue}V^{Blue}SM^{Blue}NT^{Blue}DN^{Blue}
 LM^{Blue}R^{Blue}IF^{Blue}TM^{Blue}GI^{Blue}DL^{Blue}TK^{Blue}PF^{Blue}EY^{Blue}LF^{Blue}AT^{Blue}GN^{Blue}LR^{Blue}SK^{Blue}TGL^{Blue}GLL^{Blue}Q^{Blue}DSL^{Blue}CV^{Blue}VAD^{Blue}KL^{Blue}NF^{Blue}IR^{Blue}YL^{Blue}SH^{Blue}F^{Blue}RC^{Blue}V^{Blue}HR^{Blue}GAD^{Blue}FA^{Blue}KM^{Blue}RT^{Blue}T^{Blue}
 TV^{Blue}RR^{Blue}LL^{Blue}PES^{Blue}WG^{Blue}FL^{Blue}CP^{Blue}VHT^{Blue}PD^{Blue}GE^{Blue}PC^{Blue}GL^{Blue}M^{Blue}N^{Blue}HL^{Blue}TAV^{Blue}CE^{Blue}V^{Blue}VT^{Blue}Q^{Blue}F^{Blue}V^{Blue}Y^{Blue}TA^{Blue}SIP^{Blue}ALL^{Blue}CN^{Blue}LG^{Blue}VT^{Blue}ID^{Blue}GAP^{Blue}HR^{Blue}SY^{Blue}SEC^{Blue}
 YP^{Blue}VLL^{Blue}DG^{Blue}VM^{Blue}GV^{Blue}WD^{Blue}KDL^{Blue}AP^{Blue}GI^{Blue}AD^{Blue}SL^{Blue}RH^{Blue}FK^{Blue}VL^{Blue}RE^{Blue}KR^{Blue}IP^{Blue}PW^{Blue}ME^{Blue}V^{Blue}VL^{Blue}IP^{Blue}MT^{Blue}GK^{Blue}PS^{Blue}LY^{Blue}PL^{Blue}FL^{Blue}FT^{Blue}TP^{Blue}CRL^{Blue}V^{Blue}RP^{Blue}
 VQ^{Blue}NL^{Blue}AL^{Blue}GK^{Blue}EEL^{Blue}IG^{Blue}TM^{Blue}EQ^{Blue}IF^{Blue}M^{Blue}NVA^{Blue}IF^{Blue}E^{Blue}DEV^{Blue}FAG^{Blue}VT^{Blue}TH^{Blue}QEL^{Blue}FPH^{Blue}SLL^{Blue}SV^{Blue}IAN^{Blue}FIP^{Blue}FS^{Blue}DHN^{Blue}QSP^{Blue}R^{Blue}MY^{Blue}QC^{Blue}QM^{Blue}
 GK^{Blue}Q^{Blue}TM^{Blue}GF^{Blue}PLL^{Blue}TY^{Blue}QDR^{Blue}SD^{Blue}NKL^{Blue}Y^{Blue}RL^{Blue}Q^{Blue}TP^{Blue}QS^{Blue}PL^{Blue}VR^{Blue}PS^{Blue}MY^{Blue}DY^{Blue}Y^{Blue}DM^{Blue}D^{Blue}N^{Blue}Y^{Blue}IG^{Blue}TNA^{Blue}IV^{Blue}AV^{Blue}IS^{Blue}Y^{Blue}TG^{Blue}Y^{Blue}DM^{Blue}E^{Blue}DA^{Blue}MI^{Blue}VN^{Blue}
 KAS^{Blue}WER^{Blue}GF^{Blue}AH^{Blue}GS^{Blue}VY^{Blue}KSE^{Blue}FID^{Blue}LSE^{Blue}IK^{Blue}Q^{Blue}GD^{Blue}SSL^{Blue}V^{Blue}FG^{Blue}IK^{Blue}PG^{Blue}D^{Blue}PR^{Blue}VL^{Blue}Q^{Blue}KL^{Blue}DD^{Blue}D^{Blue}GL^{Blue}PF^{Blue}IG^{Blue}AK^{Blue}L^{Blue}Q^{Blue}Y^{Blue}GD^{Blue}PY^{Blue}YS^{Blue}YL^{Blue}
 NL^{Blue}NT^{Blue}GES^{Blue}F^{Blue}VM^{Blue}Y^{Blue}Y^{Blue}SK^{Blue}ENC^{Blue}V^{Blue}VD^{Blue}NIK^{Blue}V^{Blue}CS^{Blue}ND^{Blue}TG^{Blue}SG^{Blue}K^{Blue}FC^{Blue}VC^{Blue}IT^{Blue}MR^{Blue}V^{Blue}PR^{Blue}NP^{Blue}IG^{Blue}DK^{Blue}F^{Blue}AS^{Blue}R^{Blue}H^{Blue}G^{Blue}Q^{Blue}K^{Blue}G^{Blue}IL^{Blue}S^{Blue}RL^{Blue}WP^{Blue}
 AED^{Blue}MP^{Blue}FT^{Blue}ES^{Blue}GM^{Blue}VP^{Blue}D^{Blue}IL^{Blue}FN^{Blue}PH^{Blue}GF^{Blue}PS^{Blue}RM^{Blue}TIG^{Blue}ML^{Blue}IES^{Blue}MAG^{Blue}KS^{Blue}AAL^{Blue}H^{Blue}GL^{Blue}CH^{Blue}DA^{Blue}TP^{Blue}PF^{Blue}IF^{Blue}SE^{Blue}EN^{Blue}SA^{Blue}LEY^{Blue}F^{Blue}GEM^{Blue}LK^{Blue}
 AAG^{Blue}YN^{Blue}F^{Blue}Y^{Blue}G^{Blue}TER^{Blue}L^{Blue}Y^{Blue}SG^{Blue}IS^{Blue}GLE^{Blue}E^{Blue}LE^{Blue}AD^{Blue}IF^{Blue}IG^{Blue}V^{Blue}Y^{Blue}Y^{Blue}Q^{Blue}RL^{Blue}R^{Blue}HM^{Blue}V^{Blue}SD^{Blue}K^{Blue}F^{Blue}Q^{Blue}V^{Blue}RT^{Blue}T^{Blue}GAR^{Blue}DR^{Blue}VT^{Blue}N^{Blue}Q^{Blue}PI^{Blue}GGR^{Blue}N^{Blue}V^{Blue}Q^{Blue}GG^{Blue}IR^{Blue}
 FG^{Blue}EMER^{Blue}DALL^{Blue}AH^{Blue}GT^{Blue}S^{Blue}FL^{Blue}L^{Blue}HD^{Blue}RL^{Blue}FN^{Blue}CS^{Blue}DR^{Blue}SV^{Blue}AH^{Blue}VC^{Blue}VK^{Blue}CG^{Blue}SLL^{Blue}SP^{Blue}LE^{Blue}K^{Blue}PP^{Blue}PS^{Blue}WS^{Blue}AM^{Blue}R^{Blue}N^{Blue}R^{Blue}K^{Blue}Y^{Blue}NC^{Blue}T^{Blue}LC^{Blue}SR^{Blue}SD^{Blue}
 T^{Blue}ID^{Blue}TV^{Blue}SV^{Blue}P^{Blue}Y^{Blue}V^{Blue}FR^{Blue}Y^{Blue}FVA^{Blue}E^{Blue}LA^{Blue}AM^{Blue}NIK^{Blue}V^{Blue}KL^{Blue}D^{Blue}V^{Blue}
 SG^{Green}P^{Green}TR^{Green}TR^{Green}PLE^{Green}ME^{Green}S^{Green}DE^{Green}SL^{Green}PAME^{Green}IE^{Green}CR^{Green}IT^{Green}GT^{Green}L^{Green}NG^{Green}VE^{Green}FEL^{Green}VGG^{Green}GEG^{Green}TP^{Green}EQ^{Green}GR^{Green}MT^{Green}N^{Green}K^{Green}M^{Green}ST^{Green}K^{Green}GAL^{Green}TF^{Green}SP^{Green}YLL^{Green}
 SH^{Green}VM^{Green}GY^{Green}GF^{Green}YH^{Green}FG^{Green}TY^{Green}PS^{Green}GY^{Green}EN^{Green}P^{Green}FL^{Green}HA^{Green}IN^{Green}GG^{Green}YT^{Green}N^{Green}TR^{Green}IE^{Green}KY^{Green}ED^{Green}GG^{Green}VL^{Green}H^{Green}V^{Green}S^{Green}F^{Green}SY^{Green}YE^{Green}AGR^{Green}V^{Green}IG^{Green}DF^{Green}K^{Green}VM^{Green}GT^{Green}GF^{Green}
 PED^{Green}SV^{Green}IF^{Green}TD^{Green}KI^{Green}IR^{Green}S^{Green}NAT^{Green}VE^{Green}HL^{Green}HP^{Green}MG^{Green}D^{Green}N^{Green}LD^{Green}GS^{Green}F^{Green}TR^{Green}TF^{Green}SL^{Green}RD^{Green}GG^{Green}Y^{Green}SS^{Green}V^{Green}VD^{Green}SH^{Green}M^{Green}HF^{Green}KS^{Green}AI^{Green}HP^{Green}S^{Green}IL^{Green}Q^{Green}NG^{Green}GP^{Green}
 MF^{Green}A^{Green}FR^{Green}VE^{Green}ED^{Green}HS^{Green}NTEL^{Green}G^{Green}IVE^{Green}Y^{Green}QH^{Green}AF^{Green}K^{Green}TP^{Green}DAD^{Green}AGE^{Green}ERV^{Green}


Restriction Sites: SgfI-RsrII

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN:	NM_001282779
ORF Size:	2988 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NM_001282779.2
RefSeq Size:	5356 bp
RefSeq ORF:	2991 bp
Locus ID:	84172
UniProt ID:	Q9H9Y6
Cytogenetics:	2q14.1
Protein Families:	Transcription Factors
Protein Pathways:	Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase
MW:	113.3 kDa

Gene Summary:

Eukaryotic RNA polymerase I (pol I) is responsible for the transcription of ribosomal RNA (rRNA) genes and production of rRNA, the primary component of ribosomes. Pol I is a multisubunit enzyme composed of 6 to 14 polypeptides, depending on the species. Most of the mass of the pol I complex derives from the 2 largest subunits, Rpa1 and Rpa2 in yeast. POLR1B is homologous to Rpa2 (Seither and Grummt, 1996 [PubMed 8921381]).[supplied by OMIM, Mar 2008]