

## Product datasheet for **RK300005**

### let-7c-2 Rat qPCR Template Standard (MI0000831)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGAGGTAGTAGGTTGTATGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000776
Synonyms:	let-7c-2 rno-let-7c-2; MI0000831; MIMAT0000776; rno-let-7c-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)