

## Product datasheet for **RK300006**

### let-7d Rat qPCR Template Standard (MI0000601)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** AGAGGTAGTAGGTTGCATAG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0000562

**Synonyms:** let-7d rno-let-7d; MI0000601; MIMAT0000562; rno-let-7d-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)