

Product datasheet for **RK300007**

let-7e Rat qPCR Template Standard (MI0000832)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGAGGTAGGAGGTTGTATAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000777
Synonyms:	let-7e rno-let-7e; MI0000832; MIMAT0000777; rno-let-7e-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)