

Product datasheet for **RK300024**

MIR124-3 Rat qPCR Template Standard (MI0000892)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CGTGTTACAGCGGACCTTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0004728
Synonyms:	MI0000892; MIMAT0004728; MIR124-3 rno-mir-124-3; rno-miR-124-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)