

## Product datasheet for **RK300039**

### MIR133a Rat qPCR Template Standard (MI0000906)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** AGCTGGTAAAATGGAACCAA

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0017124

**Synonyms:** MI0000906; MIMAT0017124; MIR133a rno-mir-133a; rno-miR-133a-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)